



SIDAPHARM[®]
Ophthalmic & Medical Disposable Products

SLens_{ida}

PMMA

One-Piece
Intraocular Lens

PMMA

One-Piece Intraocular Lens




— Clinical quality PMMA with UV absorber

— YAG laser compatible

— Proven Biocompatibility

ANTERIOR CHAMBER


SDANT



Optic Design	Equiconvex
Optic Size	6.00 *
Overall Length	12.50
Haptic Design	Mod. S
Refractive Index	1.489
A-Constant	115.3
AC Depth	2.95
Diopter Range	4 points of fixation -0.00 to + 30.00


POSTERIOR CHAMBER

SDHS



Optic Design	Equiconvex
Optic Size	6.00 *
Overall Length	13.00
Haptic Design	Mod. C
Refractive Index	1.489
A-Constant	118.20
AC Depth	4.80
Diopter Range	2 Dialing Holes -10.00 to + 40.00

SDALP



Optic Design	Equiconvex
Optic Size	6.00 *
Overall Length	12.50
Haptic Design	Mod. C
Refractive Index	1.489
A-Constant	118.40
AC Depth	5.00
Diopter Range	2 Dialing Holes -10.00 to + 40.00

SIDAPHARM

21, Stageiriti & 24, Em. Fili Str.
Thessaloniki, GR- 54352, Greece, European Union
T +30 2310:906660 info@sidapharm.gr
F: +30 2310 989846 www.sidapharm.gr

