


PRODUCT SPECIFICATIONS

for **SIDA-Fluoro**

Fluorescein Sodium Ophthalmic strips

Ref. No.	91000
Intended Use	Fluorescein Sodium Ophthalmic Strip is indicated for use as an aid for contact lens fitment.
Certification	CE 0653
Classification	Class Is, according to Annex IX of Council Directive 93/42/EEC
Sterilization Method	ETO
Conformity	MEDDEV 2.4/1 Rev. 9 - June 2010 Rule 5
Applicable Standards	ISO 9001, ISO 13485, ISO 14971, ISO 14155, ISO 10993, ISO 11607, ISO 14644, ISO 11737, ISO 11140, ISO 11138, ISO 13408
Depiction	
Shelf life	5 years
Primary Packing	Pack of 100 pieces
Secondary Packing	Paper Box