

INTRODUCING THE NEW
ESC² OPHTHALMOLOGY
PLATFORM



DESIGNED TO BOOST
PRODUCTIVITY

STREAMLINE PATIENT HANDLING
IMPROVE PATIENT FLOW AND THROUGHPUT
ENHANCE PATIENT, SURGEON AND STAFF SATISFACTION

ESC² DESIGNED BY CLINICIANS FOR CLINICIANS

THE MOST ADVANCED EYE CHAIR IN TODAY'S MARKET



Patented Articulating Head Rest with Infinite Positioning

Contoured, Multi-density Pad with DURA-GEL[®] Insert for Greater Patient Comfort

Integrated Powered Leg Extension for Taller Patients

Superior Surgeon Leg Clearance

Low Profile Base for Smaller Spaces and Greater Clearance

SMART CASTER TECHNOLOGY with Auto Braking and Advanced Steering for Improved Patient and Staff Safety

ESC²

ADVANCED SINGLE PATIENT PLATFORM

600 POUNDS
600 lb (272 kg)
Weight Capacity
EVENLY DISTRIBUTED
Highest on the
Market Today

ESC² BUILT TODAY. DESIGNED FOR TOMORROW



Positioning at Your Fingertips:
Intuitive Hand Pendant with Pre-Sets
and 3 Memory Positions

Integrated Ergonomic
Push Handles for
Increased Staff Safety



Wider Patient Surface
- 26 inch (66 cm) for
Enhanced Patient
Comfort and Safety

Counter Rotating Side Rails
with Urethane Comfort Arm
Rests for Improved Patient
Safety and Comfort



Two-position Foot Rest For
Patient Comfort and Safety

Industry Leading Height Range –
20 to 40 inches (51 to 102 cm) for
Greater Procedure Flexibility

ESC²
THE NEXT GENERATION
PLATFORM

ESC² IMPROVED VERSATILITY, EFFICIENCY, AND COMFORT



Optional Accessory:
Full "U" Wrist Rest
Item No: HSA078600



0-15 degrees
Seat Tilt
for Patient
Comfort and
Safety

Optional Accessory:
O₂ Tank Holder
Item No: HSA080009

Complete the ESC²
with the new Hausted
Surgical Stool
Item No: HSS500



Specifications:

	ESC ² EYEST	
	inch	cm
Width:	32	81
Length:	83	211
Height: maximum	40	102
Height: minimum	20	51
Patient Surface Width:	26	66
Patient Surface Length: minimum	75	191
Patient Surface Length: maximum	83	211
Side Rail Length:	18.5	47
Side Rail Height:	10.25	26
Caster Diameter:	5	13
Multi-Density Sectional Pad: Dura-Gel™ insert, urethane cover	3	8
	degree °1	
Trendelenburg:	15	
Seat Tilt:	0 - 15	
Backrest Range:	0 - 85	
	volts	
Dual Power:	110 - 240	
	lb	kg
Product Weight:	265	120
Weight Capacity: <small>EVENLY DISTRIBUTED</small>	600	272
Shipping Weight:	315	143

Standard Accessories:

Description	UOM
Telescoping IV Pole	1/ea
Safety Strap	1/ea
Utility Bag Hook	1/ea

Optional Accessories:

Item No.	Description	UOM
HSA078600	Full "U" Wrist Rest	1/ea
HSA078500	Dual Wrist Rest	1/ea
HSA080003	OxyFlex II Oxygen Delivery System	1/ea

Item No.	Description	UOM
HSA008000	OxyFlex II Diffusion Trays	25/cs
H0103ESC2	Contoured Head Rests	3/set
HSA080014	Hi/Lo Foot Switch	1/ea

Item No.	Description	UOM
HSA080013	Shoulder Support Rails	1/pr
HSA080020	Foot End Push Handles	1/pr
HSA080009	O ₂ Tank Holder	1/ea

WARNING: Cancer and Reproductive Harm -
www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



CE Partner4U BV
Esdoolaan 13
3951 DB Maarn
The Netherlands
www.cephartner4u.eu



© 2021, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



GF2000176/Rev/E21