


**Sterimedix**<sup>®</sup>  
ophthalmic

# EXPERTISE

in single-use ophthalmic surgical products




[www.sterimedix.com](http://www.sterimedix.com)



# Consistency in a Changing World

manufacturing single-use ophthalmic cannulas and handpieces

With over 30 years' experience as a manufacturer of single use ophthalmic surgical products, Sterimedix is now the leading European manufacturer of single use ophthalmic cannulas.




meeting the demands of modern ophthalmic surgery techniques



## Index

### PAGE

|    |                                                          |
|----|----------------------------------------------------------|
| 4  | Refractive Cannulas                                      |
| 4  | Scleral Markers                                          |
| 5  | Anaesthesia Needles & Cannulas                           |
| 6  | Anterior Chamber Cannulas & Maintainers                  |
| 6  | Transzonular Injection Cannula                           |
| 7  | Viscoelastic Cannulas                                    |
| 8  | Cystotomes                                               |
| 10 | Hydrodissection Cannulas                                 |
| 11 | Hydrodelineation Cannulas                                |
| 11 | Cortical Cleaving Cannulas                               |
| 11 | Viscoexpression Cannulas                                 |
| 12 | Nucleus Removal Cannulas                                 |
| 12 | Cortex Removal Cannulas                                  |
| 13 | Capsule Polishing Cannulas                               |
| 14 | Irrigating / Aspirating Cannulas and Handpieces          |
| 15 | Simcoe Irrigating / Aspirating Cannulas                  |
| 15 | Irrigating Microsurgical Handle                          |
| 16 | Bimanual Irrigating / Aspirating Cannulas and Handpieces |
| 17 | Bimanual I/A Systems                                     |
| 18 | Lachrymal Cannulas                                       |
| 19 | Trocar Sets                                              |
| 20 | Vitreoretinal Cannulas                                   |
| 21 | Vitreoretinal Soft Tips & Brushes                        |
| 22 | Backflush Handles & Combination Sets                     |
| 23 | Heavy Liquid Infusion Handles                            |
| 23 | Charles Flute Handles & Combination Sets                 |
| 23 | Information for Kit Packers                              |
| 24 | Product Index                                            |



Products are supplied in different formats to meet the needs of our customers:

- ▶ Packed sterile in convenient box quantities of either 5, 10 or 100.
- ▶ Packed sterile in bulk for our industry partners
- ▶ Packed non-sterile, in bulk for kit packers.  
Most non-sterile products are packed in a protective sheath labelled with the product code & LOT number.



### QUALITY ASSURANCE


Sterimedix Limited is committed to producing medical devices to the highest possible quality. The Company's Quality Management System is accredited to ISO 13485:2016. Our products are CE marked, registered with the United States FDA and have appropriate approvals in all other markets where Sterimedix is active.

# Refractive Cannulas / Scleral Markers

## Refractive Cannulas

For irrigating either the anterior or posterior sides of the flap and the stroma in laser refractive procedures.

|       |                                                                                          |                             |                                                                |                                                        |                                                                                       | BOX |     |
|-------|------------------------------------------------------------------------------------------|-----------------------------|----------------------------------------------------------------|--------------------------------------------------------|---------------------------------------------------------------------------------------|-----|-----|
|       |                                                                                          |                             |                                                                |                                                        |                                                                                       |     |     |
| M2804 | 25 G    | 25 G x 7/8"<br>0.5 x 22mm   | Irrigating lasik cannula.<br>Formed 7mm from tip               | 6 x 0.15mm ports<br>at 2, 4 and 6mm<br>from tip        |    | 10  | 100 |
| M2817 | 27 G    | 27 G x 7/8"<br>0.4 x 22mm   | Irrigating lasik cannula.<br>Formed 9mm from tip               | 4 x 0.15mm ports<br>at 2 and 4mm<br>from tip           |    | 10  | 100 |
| M2822 | 25 G    | 25 G x 7/8"<br>0.5 x 22mm   | Irrigating lasik cannula.<br>Formed 7mm from tip               | 6 x 0.15mm ports<br>at 2, 4 and 6mm<br>from closed tip |    | 10  | 100 |
| M2824 | 26 G    | 26 G x 7/8"<br>0.45 x 22mm  | Irrigating lasik cannula.<br>Formed 9.5mm from tip.<br>(Slade) | Spatulated tip                                         |    | 10  | 100 |
| M2801 | 25 G   | 25 G x 7/8"<br>0.5 x 22mm   | Irrigating lasik cannula.<br>Angled. (Lin)                     | Flattened tip<br>2 x 0.25mm ports                      |   | 10  | 100 |
| M2802 | 25 G  | 25 G x 1 1/8"<br>0.5 x 29mm | Irrigating lasik cannula.<br>Formed. (Buratto)                 | 2 x 0.2mm ports                                        |  | 10  | 100 |
| M3914 | 23 G  | 23 G x 7/8"<br>0.64 x 22mm  | Irrigating lasik cannula.<br>Formed 10.5mm from tip            | Flattened tip                                          |  | 10  | 100 |



## Scleral Markers (Sterile and Non Sterile)

Used to mark the point for any pars plana incision or injection point

| M7500 |  | 3.5mm / 4.0mm | Scleral marker. 80mm long         |                                                                                       | BOX |
|-------|-------------------------------------------------------------------------------------|---------------|-----------------------------------|---------------------------------------------------------------------------------------|-----|
|       |                                                                                     |               |                                   |                                                                                       | 10  |
| M7510 |  | 3.5mm         | Scleral fixation plate. 75mm long |  | 10  |

# Anaesthetic Cannulas and Needles

## Anaesthetic Cannulas

For administering anaesthetic into the sub tenon space.

| M0278 | 19 G |   | 19 G x 1"<br>1.1 x 25mm | Sub tenon cannula. Curved | Flattened tip |    | BOX |     |
|-------|------|-------------------------------------------------------------------------------------|-------------------------|---------------------------|---------------|---------------------------------------------------------------------------------------|-----|-----|
|       |      |                                                                                     |                         |                           |               |                                                                                       | 10  | 100 |
| M0279 | 20 G |  | 20 G x 1"<br>0.9 x 25mm | Sub tenon cannula. Curved | Flattened tip |  | BOX |     |
|       |      |                                                                                     |                         |                           |               |                                                                                       | 10  | 100 |

## Anaesthetic Needles

For administering anaesthetic around the globe.

| M0275  | 25 G |  | 25 G x 1 1/2"<br>0.5 x 38mm  | Retrobulbar. Tri-facet |  | BOX |     |
|--------|------|-------------------------------------------------------------------------------------|------------------------------|------------------------|---------------------------------------------------------------------------------------|-----|-----|
|        |      |                                                                                     |                              |                        |                                                                                       | 10  | 100 |
| M0637  | 25 G |  | 25 G x 1 1/2"<br>0.5 x 38mm  | Retrobulbar (Atkinson) |  | BOX |     |
|        |      |                                                                                     |                              |                        |                                                                                       | 10  | 100 |
| M0638  | 23 G |  | 23 G x 1 1/2"<br>0.64 x 38mm | Retrobulbar (Atkinson) |  | BOX |     |
|        |      |                                                                                     |                              |                        |                                                                                       | 10  | 100 |
| M0641A | 27 G |  | 27 G x 1 1/4"<br>0.4 x 32mm  | Peribulbar (Atkinson)  |  | BOX |     |
|        |      |                                                                                     |                              |                        |                                                                                       | 10  | 100 |
| M0642  | 25 G |  | 25 G x 1 1/4"<br>0.5 x 32mm  | Peribulbar (Atkinson)  |  | BOX |     |
|        |      |                                                                                     |                              |                        |                                                                                       | 10  | 100 |
| M0643  | 23 G |  | 23 G x 1 1/4"<br>0.64 x 32mm | Peribulbar (Atkinson)  |  | BOX |     |
|        |      |                                                                                     |                              |                        |                                                                                       | 10  | 100 |
| M0641  | 27 G |  | 27 G x 7/8"<br>0.4 x 22mm    | Peribulbar (Atkinson)  |  | BOX |     |
|        |      |                                                                                     |                              |                        |                                                                                       | 10  | 100 |
| M0642A | 25 G |  | 25 G x 7/8"<br>0.5 x 22mm    | Peribulbar (Atkinson)  |  | BOX |     |
|        |      |                                                                                     |                              |                        |                                                                                       | 10  | 100 |

# Anterior Chamber Maintenance

## Anterior Chamber Maintainers

Self retaining cannulas for maintaining and forming the anterior or posterior chamber by infusing or aspirating balanced salt solution (BSS), air, or viscoelastic.


|       |      |                                                                                   |                              |                                                                                                                      |                                                                                     |     |
|-------|------|-----------------------------------------------------------------------------------|------------------------------|----------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----|
| M2658 | 23 G |  | 23 G x 3/64"<br>0.64 x 3.5mm | Anterior chamber maintainer.<br>150mm silicone tubing.<br>Annular rings allow cannula to be self retaining (Lewicky) |  | BOX |
|       |      |                                                                                   |                              |                                                                                                                      |                                                                                     | 10  |
| M2657 | 20 G |  | 20 G x 3/64"<br>0.64 x 3.5mm | Anterior chamber maintainer.<br>150mm silicone tubing.<br>Annular rings allow cannula to be self retaining (Lewicky) |  | 10  |
|       |      |                                                                                   |                              |                                                                                                                      |                                                                                     |     |

## Anterior Chamber Maintenance Cannulas

For maintaining and forming the anterior or posterior chamber by infusing or aspirating balanced salt solution (BSS), air, or viscoelastic.

| M2280  | 27 G |    | 27 G x 7/8"<br>0.4 x 22mm  | Transzonular injection cannula. Angled      | Angled 33°, 2mm from tip                                              |    | BOX |     |
|--------|------|-------------------------------------------------------------------------------------|----------------------------|---------------------------------------------|-----------------------------------------------------------------------|---------------------------------------------------------------------------------------|-----|-----|
|        |      |                                                                                     |                            |                                             |                                                                       |                                                                                       | 10  | 100 |
| M2639  | 25 G |   | 25 G x 1"<br>0.5 x 25mm    | Anterior chamber cannula. Angled. (Rycroft) | Angled 35°, 7.5mm from tip                                            |   | 10  | 100 |
| M2270  | 30 G |  | 30 G x 7/8"<br>0.3 x 22mm  | Anterior chamber cannula. Angled. (Rycroft) | Angled 35°, 4mm from tip                                              |  | 10  | 100 |
| M2270A | 30 G |  | 30 G x 7/8"<br>0.3 x 22mm  | Anterior chamber cannula. Angled (Rycroft)  | Angled 35°, 8mm from tip                                              |  | 10  | 100 |
| M2273  | 30 G |  | 30 G x 7/8"<br>0.3 x 22mm  | Anterior chamber cannula. Angled. (Rycroft) | Angled 35°, 5.5mm from tip                                            |  | 10  | 100 |
| M2273A | 27 G |  | 27 G x 7/8"<br>0.4 x 22mm  | Anterior chamber cannula. Angled. (Rycroft) | Angled 35°, 4.5mm from tip                                            |  | 10  | 100 |
| M2279  | 27 G |  | 27 G x 7/8"<br>0.4 x 22mm  | Anterior chamber cannula. Angled. (Rycroft) | Angled 45°, 11mm from tip                                             |  | 10  | 100 |
| M2275  | 25 G |  | 25 G x 7/8"<br>0.5 x 22mm  | Anterior chamber cannula. Angled. (Rycroft) | Angled 35°, 5mm from tip                                              |  | 10  | 100 |
| M2273G | 23 G |  | 23 G x 7/8"<br>0.64 x 22mm | Anterior chamber cannula. Angled. (Rycroft) | Angled 35°, 5.5mm from tip                                            |  | 10  | 100 |
| M2642  | 21 G |  | 21 G x 7/8"<br>0.81 x 22mm | Anterior chamber cannula. Angled. (Rycroft) | Angled 35°, 7.5mm from tip                                            |  | 10  | 100 |
| M2273B | 20 G |  | 20 G x 7/8"<br>0.9 x 22mm  | Anterior chamber cannula. Angled. (Rycroft) | Angled 45°, 7.5mm from tip                                            |  | 10  | 100 |
| M9280  | 19 G |  | 19 G x 1"<br>1.1 x 25mm    | Washout cannula. Angled. (Bishop Harmon)    | Angled 32°, 8mm from tip. Tip flattened in the horizontal plane       |  | 10  | 100 |
| M9281  | 16 G |  | 16 G x 1"<br>1.6 x 25mm    | Washout cannula. Angled                     | Angled 30°, 6.5mm from the tip. Tip flattened in the horizontal plane |  | 10  | 100 |

# Viscoelastic Cannulas



## Viscoelastic Cannulas

For maintaining and forming the anterior or posterior chamber by infusing or aspirating viscoelastic.

|        |      |  |                            |                                            |                                        | BOX |        |
|--------|------|--|----------------------------|--------------------------------------------|----------------------------------------|-----|--------|
|        |      |  |                            |                                            |                                        |     |        |
| M3926A | 27 G |  | 27 G x 7/8"<br>0.4 x 22mm  | Viscoelastic cannula.<br>Angled            | Angled 35°,<br>9mm from tip            |     | 10 100 |
| M3927  | 27 G |  | 27 G x 7/8"<br>0.4 x 22mm  | Viscoelastic cannula.<br>Angled            | Angled 25°,<br>8mm from tip            |     | 10 100 |
| M2272  | 26 G |  | 26 G x 7/8"<br>0.45 x 22mm | Viscoelastic cannula.<br>Angled            | Angled 35°,<br>9mm from tip            |     | 10 100 |
| M2639A | 25 G |  | 25 G x 7/8"<br>0.5 x 22mm  | Viscoelastic cannula.<br>Angled            | Angled 35°,<br>9mm from tip            |     | 10 100 |
| M3909  | 23 G |  | 23 G x 7/8"<br>0.64 x 22mm | Viscoelastic cannula.<br>Angled            | Angled 45°,<br>9mm from tip            |     | 10 100 |
| M3910A | 21 G |  | 21 G x 7/8"<br>0.81 x 22mm | Viscoelastic cannula.<br>Angled            | Angled 45°,<br>9mm from tip            |     | 10 100 |
| M3910  | 21 G |  | 21 G x 7/8"<br>0.81 x 22mm | Viscoelastic cannula.                      |                                        |     | 10 100 |
| M3917  | 27 G |  | 27 G x 7/8"<br>0.4 x 22mm  | Viscoelastic cannula.<br>Angled and formed | Angled<br>and formed,<br>11mm from tip |     | 10 100 |

# Cystotomes


For capsularhexis, intercapsular, endocapsular, or can opening capsulotomies during cataract surgery.

## Formed Irrigating Cystotomes

|        |                                                                                          |                           |                                |                |                                                                                       | BOX |     |
|--------|------------------------------------------------------------------------------------------|---------------------------|--------------------------------|----------------|---------------------------------------------------------------------------------------|-----|-----|
|        |                                                                                          |                           |                                |                |                                                                                       |     |     |
| M3630  | 30 G    | 30 G x 5/8"<br>0.3 x 16mm | Cystotome. Formed              |                |    | 10  | 100 |
| M3610B | 27 G    | 27 G x 5/8"<br>0.4 x 16mm | Cystotome. Formed              |                |    | 10  | 100 |
| M3610  | 25 G    | 25 G x 5/8"<br>0.5 x 16mm | Cystotome. Formed              |                |    | 10  | 100 |
| M3610S | 25 G    | 25 G x 5/8"<br>0.5 x 16mm | Cystotome. Formed              | Short tip      |    | 10  | 100 |
| M3602  | 25 G    | 25 G x 5/8"<br>0.5 x 16mm | Cystotome. Formed              | Side cutting   |    | 10  | 100 |
| M3610C | 25 G   | 25 G x 5/8"<br>0.5 x 16mm | Cystotome. Formed              | Tip angled 35° |   | 10  | 100 |
| M3627A | 27 G  | 27 G x 5/8"<br>0.4 x 16mm | Cystotome. Short curve. Formed |                |  | 10  | 100 |
| M3627  | 25 G  | 25 G x 5/8"<br>0.5 x 16mm | Cystotome. Short curve. Formed |                |  | 10  | 100 |
| M3626  | 27 G  | 27 G x 5/8"<br>0.4 x 16mm | Cystotome. Angled              |                |  | 10  | 100 |
| M3626A | 25 G  | 25 G x 5/8"<br>0.5 x 16mm | Cystotome. Angled              |                |  | 10  | 100 |
| M3629  | 25 G  | 25 G x 5/8"<br>0.5 x 16mm | Cystotome. Formed              | Reverse tip    |  | 10  | 100 |
| M3615  | 27 G  | 27 G x 7/8"<br>0.4 x 22mm | Long cystotome. Formed         |                |  | 10  | 100 |

# Cystotomes

For capsularhexis, intercapsular, endocapsular, or can opening capsulotomies during cataract surgery.



## Formed Capsularhexis Cystotomes

|        |                                                                                          |                           |                                |                                                                                       | BOX |     |
|--------|------------------------------------------------------------------------------------------|---------------------------|--------------------------------|---------------------------------------------------------------------------------------|-----|-----|
|        |                                                                                          |                           |                                |                                                                                       |     |     |
| M3000F | 27 G  | 27 G x 5/8"<br>0.4 x 16mm | Cystotome. Formed              |  | 10  | 100 |
| M3000G | 25 G  | 25 G x 5/8"<br>0.5 x 16mm | Cystotome. Formed              |  | 10  | 100 |
| M3000E | 27 G  | 27 G x 5/8"<br>0.4 x 16mm | Cystotome. Short curve. Formed |  | 10  | 100 |
| M3000D | 25 G  | 25 G x 5/8"<br>0.5 x 16mm | Cystotome. Short curve. Formed |  | 10  | 100 |

## Formed Endocapsular Cystotome

|                            |                                                                                          |                           |                                         |                                                                                       | BOX |     |
|----------------------------|------------------------------------------------------------------------------------------|---------------------------|-----------------------------------------|---------------------------------------------------------------------------------------|-----|-----|
|                            |                                                                                          |                           |                                         |                                                                                       |     |     |
| <b>NEW</b><br><b>M3001</b> | 27 G  | 27 G x 5/8"<br>0.4 x 16mm | Endocapsular cystotome. Formed. (Joyce) |  | 10  | 100 |

# Hydrodissection Cannulas

For delivering fluids to facilitate the separation of the cortex from the nucleus & capsule.

## Hydrodissection Cannulas

| Model  | Gauge | Image | Dimensions                | Design                                | Tip Description                                                           | Image | BOX |     |
|--------|-------|-------|---------------------------|---------------------------------------|---------------------------------------------------------------------------|-------|-----|-----|
|        |       |       |                           |                                       |                                                                           |       |     |     |
| M3919  | 27 G  |       | 27 G x 7/8"<br>0.4 x 22mm | Hydrodissector.<br>Angled. (Akahoshi) | 2mm tapered tip                                                           |       | 10  | 100 |
| M2273C | 27 G  |       | 27 G x 7/8"<br>0.4 x 22mm | Hydrodissector.<br>Curved. (Sauter)   | Tip flattened in vertical plane                                           |       | 10  | 100 |
| M2273D | 25 G  |       | 25 G x 7/8"<br>0.5 x 22mm | Hydrodissector.<br>Curved. (Sauter)   | Tip flattened in vertical plane                                           |       | 10  | 100 |
| M2273E | 27 G  |       | 27 G x 7/8"<br>0.4 x 22mm | Hydrodissector.<br>Curved. (Helsinki) | Tip flattened in horizontal plane                                         |       | 10  | 100 |
| M2273F | 25 G  |       | 25 G x 7/8"<br>0.5 x 22mm | Hydrodissector.<br>Curved. (Helsinki) | Tip flattened in horizontal plane                                         |       | 10  | 100 |
| M2274  | 25 G  |       | 25 G x 7/8"<br>0.5 x 22mm | Hydrodissector.<br>Curved. (Kellan)   | Tip flattened in vertical plane                                           |       | 10  | 100 |
| M3992  | 30 G  |       | 30 G x 7/8"<br>0.3 x 22mm | Hydrodissector. Angled                | Angled 7mm from tip.<br>Tip tapered and flattened in the horizontal plane |       | 10  | 100 |
| M3902  | 27 G  |       | 27 G x 7/8"<br>0.4 x 22mm | Hydrodissector. Angled                | Angled 9mm from tip.<br>Tip tapered and flattened in the horizontal plane |       | 10  | 100 |
| M3900  | 25 G  |       | 25 G x 7/8"<br>0.5 x 22mm | Hydrodissector. Angled                | Angled 7mm from tip.<br>Tip tapered and flattened in the horizontal plane |       | 10  | 100 |
| M2273H | 27 G  |       | 27 G x 7/8"<br>0.4 x 22mm | Hydrodissector. Curved                | Rounded tip                                                               |       | 10  | 100 |

# Hydrodissection / Hydrodelineation / Viscoexpression Cannulas

For delivering fluids to facilitate the separation of the cortex from the nucleus & capsule.

## Hydrodissection Cannulas

| Model     | Gauge | Image                                                                             | Size                       | Type                                | Description                                                                        | Image                                                                               | BOX |     |
|-----------|-------|-----------------------------------------------------------------------------------|----------------------------|-------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----|-----|
|           |       |                                                                                   |                            |                                     |                                                                                    |                                                                                     | 10  | 100 |
| M3904     | 27 G  |  | 27 G x 3/8"<br>0.4 x 19mm  | Hydrodissector. (Chang)             | Tip angled 90° and flattened in the vertical plane                                 |  | 10  | 100 |
| M3912     | 23 G  |  | 23 G x 7/8"<br>0.64 x 22mm | Hydrodissector. Angled. (Mackool)   | Angled 30° 8mm from tip. Flattened 4mm from tip. Tip flattened in horizontal plane |  | 10  | 100 |
| M3906A    | 25 G  |  | 25 G x 7/8"<br>0.5 x 22mm  | Hydrodissector                      | Round tip, 1.75mm wide. J-shaped tip                                               |  | 10  | 100 |
| NEW M3920 | 27 G  |  | 27 G x 7/8"<br>0.4 x 22mm  | Hydrodissector. Curved. (Frankfurt) |                                                                                    |  | 10  | 100 |
| NEW M3922 | 27 G  |  | 27 G x 7/8"<br>0.4 x 22mm  | Hydrodissector. Formed. (Daya)      | Tip flattened in horizontal plane                                                  |  | 10  | 100 |

## Cortical Cleaving Cannulas

| Model  | Gauge | Image                                                                               | Size                      | Type                                     | Description                                                                 | Image                                                                                 | BOX |     |
|--------|-------|-------------------------------------------------------------------------------------|---------------------------|------------------------------------------|-----------------------------------------------------------------------------|---------------------------------------------------------------------------------------|-----|-----|
|        |       |                                                                                     |                           |                                          |                                                                             |                                                                                       | 10  | 100 |
| M3908  | 25 G  |  | 25 G x 7/8"<br>0.5 x 22mm | Cortical cleaving hydrodissector. Angled | Tapered to 0.3mm at tip. Angled at 45°, 1mm from tip                        |  | 10  | 100 |
| M3911  | 27 G  |  | 27 G x 7/8"<br>0.4 x 22mm | Cortical cleaving hydrodissector. Angled | Angled at 45°, 11mm from tip. Flattened in horizontal plane 2.75mm from tip |  | 10  | 100 |
| M3911A | 27 G  |  | 27 G x 7/8"<br>0.4 x 22mm | Cortical cleaving hydrodissector. Angled | Angled 45°, 11mm from round tip                                             |  | 10  | 100 |

## Hydrodelineation Cannulas

| Model     | Gauge | Image                                                                               | Size                      | Type                    | Description                                            | Image                                                                                 | BOX |     |
|-----------|-------|-------------------------------------------------------------------------------------|---------------------------|-------------------------|--------------------------------------------------------|---------------------------------------------------------------------------------------|-----|-----|
|           |       |                                                                                     |                           |                         |                                                        |                                                                                       | 10  | 100 |
| M3907     | 25 G  |  | 25 G x 7/8"<br>0.5 x 22mm | Hydrodelineator         | Tapered to 0.3mm at tip. Angled at 45°, 4.5mm from tip |  | 10  | 100 |
| NEW M3921 | 25 G  |  | 25 G x 7/8"<br>0.5 x 22mm | Hydrodelineator. Curved | Tapered to 0.3mm at tip                                |  | 10  | 100 |

## Viscoexpression Cannulas

| Model  | Gauge | Image                                                                               | Size                      | Type                    | Description                               | Image                                                                                 | BOX |     |
|--------|-------|-------------------------------------------------------------------------------------|---------------------------|-------------------------|-------------------------------------------|---------------------------------------------------------------------------------------|-----|-----|
|        |       |                                                                                     |                           |                         |                                           |                                                                                       | 10  | 100 |
| M3899  | 25 G  |  | 25 G x 7/8"<br>0.5 x 22mm | Viscoexpression cannula | Flat tip angled 2mm from tip (Corydon)    |  | 10  | 100 |
| M3899A | 25 G  |  | 25 G x 7/8"<br>0.5 x 22mm | Viscoexpression cannula | Round tip angled 1.8mm from tip (Corydon) |  | 10  | 100 |

# Nucleus / Cortex Removal Cannulas

## Nucleus Removal Cannulas

For the aspiration of cortex

|              |                                                                                               |                             |                   |                          |                                                                                     |     |     |
|--------------|-----------------------------------------------------------------------------------------------|-----------------------------|-------------------|--------------------------|-------------------------------------------------------------------------------------|-----|-----|
| <b>M4619</b> | <b>25 G</b>  | 25 G x 1 1/2"<br>0.5 x 38mm | Irrigating vectis | 5mm wide.<br>Single port |  | BOX |     |
|              |                                                                                               |                             |                   |                          |                                                                                     | 10  | 100 |
| <b>M4620</b> | <b>25 G</b>  | 25 G x 1 1/2"<br>0.5 x 38mm | Irrigating vectis | 5mm wide.<br>Twin ports  |  | 10  | 100 |

## Cortex Removal Cannulas

For either removing the nucleus, large fragments of the nucleus, large pieces of cortical debris, or aspirating cortical debris.

|                       |                                                                                                 |                              |                                         |                            |                                                                                       |     |     |
|-----------------------|-------------------------------------------------------------------------------------------------|------------------------------|-----------------------------------------|----------------------------|---------------------------------------------------------------------------------------|-----|-----|
| <b>M6576</b>          | <b>23 G</b>    | 23 G x 7/8"<br>0.64 x 22mm   | Cortex extractor. Curved                | 0.3mm top port             |    | BOX |     |
|                       |                                                                                                 |                              |                                         |                            |                                                                                       | 10  | 100 |
| <b>M6575</b>          | <b>23 G</b>    | 23 G x 5/8"<br>0.64 x 16mm   | Cortex extractor. Curved                | 0.3mm top port             |    | 10  | 100 |
| <b>M6575A</b>         | <b>23 G</b>  | 23 G x 5/8"<br>0.64 x 16mm   | Cortex extractor. Curved                | 0.4mm top port             |  | 10  | 100 |
| <b>NEW<br/>M6575B</b> | <b>23 G</b>  | 23 G x 5/8"<br>0.64 x 16mm   | Cortex extractor. Curved                | 0.3mm top port             |  | 10  | 100 |
| <b>M6576A</b>         | <b>23 G</b>  | 23 G x 5/8"<br>0.64 x 16mm   | Cortex extractor. Curved                | 0.3mm top port             |  | 10  | 100 |
| <b>M6594</b>          | <b>23 G</b>  | 23 G x 1/2"<br>0.64 x 12.5mm | Cortex extractor.<br>Curved. (Charleux) | 0.3mm bevelled end<br>port |  | 10  | 100 |
| <b>M6596</b>          | <b>24 G</b>  | 24 G x 1/2"<br>0.56 x 12.5mm | Cortex extractor.<br>Curved. (Charleux) | 0.3mm bevelled end<br>port |  | 10  | 100 |
| <b>M6595</b>          | <b>23 G</b>  | 23 G x 1/2"<br>0.64 x 12.5mm | Cortex extractor.<br>Curved. (Charleux) | 0.3mm bevelled end<br>port |  | 10  | 100 |

# Capsule Polishing Cannulas

## Capsule Polishing Cannulas

For polishing and cleaning the residual cortex from the capsular bag. Silicone tipped polishers provide for gentle polishing and atraumatic insertion and removal. Micro etched tips are slightly roughened to facilitate gentle polishing.

|        |                                                                                        |                            |                          |                                                                                  |                                                                                     | BOX |     |
|--------|----------------------------------------------------------------------------------------|----------------------------|--------------------------|----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----|-----|
|        |                                                                                        |                            |                          |                                                                                  |                                                                                     |     |     |
| M8604A | 23 G  | 23 G x 7/8"<br>0.64 x 22mm | Capsule polisher. Angled | Angled 9mm from tip.<br>Aspiration through 0.3mm side port.<br>Micro-etched tip  |  | 10  | 100 |
| M8604  | 21 G  | 21 G x 1"<br>0.81 x 25mm   | Capsule polisher. Angled | Angled 12.5mm from tip. Aspiration through 0.3mm side port.<br>Micro-etched tip  |  | 10  | 100 |
| M8604B | 21 G  | 21 G x 7/8"<br>0.81 x 22mm | Capsule polisher. Angled | Angled 9mm from tip.<br>Aspiration through 0.3mm side port.<br>Micro-etched tip. |  | 10  | 100 |
| M8645  | 23 G  | 23 G x 7/8"<br>0.64 x 22mm | Capsule scraper. Angled  | Angled 0.5 & 7mm from tip. Micro-etched distal end.                              |  | 10  | 100 |

## Silicone Tipped Polishers

|        |                                                                                          |                            |                                  |  |                                                                                       | BOX |     |
|--------|------------------------------------------------------------------------------------------|----------------------------|----------------------------------|--|---------------------------------------------------------------------------------------|-----|-----|
|        |                                                                                          |                            |                                  |  |                                                                                       |     |     |
| M8605  | 27 G  | 27 G x 7/8"<br>0.4 x 22mm  | Silicone tipped polisher. Angled |  |  | 10  | 100 |
| M8605A | 23 G  | 23 G x 7/8"<br>0.64 x 22mm | Silicone tipped polisher. Angled |  |  | 10  | 100 |

## Olive Tipped Cannulas

|        |                                                                                          |                          |                               |                                            |                                                                                       | BOX |     |
|--------|------------------------------------------------------------------------------------------|--------------------------|-------------------------------|--------------------------------------------|---------------------------------------------------------------------------------------|-----|-----|
|        |                                                                                          |                          |                               |                                            |                                                                                       |     |     |
| M8622  | 25 G  | 25 G x 1"<br>0.5 x 25mm  | Olive tipped polisher. Angled | Angled 10mm from tip.<br>Micro-etched tip. |  | 10  | 100 |
| M8620C | 23 G  | 23 G x 1"<br>0.64 x 25mm | Olive tipped polisher. Curved |                                            |  | 10  | 100 |
| M8621  | 23 G  | 23 G x 1"<br>0.64 x 25mm | Olive tipped polisher. Curved | Micro-etched tip                           |  | 10  | 100 |


# Irrigating/Aspirating Cannulas and Handpieces

For the aspiration of cortex while simultaneously maintaining the Anterior Chamber with continuous irrigating flow.

## Coaxial Cannulas

| M6611 | 21 G  | 21 G x 1½"<br>0.81 x 38mm | I/A Straight Cannula |  | BOX |     |
|-------|----------------------------------------------------------------------------------------|---------------------------|----------------------|-------------------------------------------------------------------------------------|-----|-----|
|       |                                                                                        |                           |                      |                                                                                     | 10  | 100 |

## Coaxial Handpieces

| M5592 | 17/24G    | 1.47 / 0.56mm | I/A handpiece for 2.2mm incisions                  | 45° angled tip |    | BOX |  |
|-------|--------------------------------------------------------------------------------------------|---------------|----------------------------------------------------|----------------|---------------------------------------------------------------------------------------|-----|--|
|       |                                                                                            |               |                                                    |                |                                                                                       | 5   |  |
| M5593 | 17/24G    | 1.47 / 0.56mm | I/A handpiece for 2.2mm incisions                  | 60° angled tip |    | 5   |  |
| M5594 | 17/24G    | 1.47 / 0.56mm | I/A handpiece for 2.2mm incisions                  | 90° angled tip |    | 5   |  |
| M5591 | 17/24G    | 1.47 / 0.56mm | I/A handpiece for 2.2mm incisions                  | Curved tip     |    | 5   |  |
| M5580 | 16/21G   | 1.6 / 0.81mm  | I/A handpiece for 2.7 / 2.8mm incisions            | Straight tip   |   | 5   |  |
| M5582 | 16/21G  | 1.6 / 0.81mm  | I/A handpiece for 2.7 / 2.8mm incisions            | 45° angled tip |  | 5   |  |
| M5583 | 16/21G  | 1.6 / 0.81mm  | I/A handpiece for 2.7 / 2.8mm incisions            | 60° angled tip |  | 5   |  |
| M5584 | 16/21G  | 1.6 / 0.81mm  | I/A handpiece for 2.7 / 2.8mm incisions            | 90° angled tip |  | 5   |  |
| M5581 | 16/21G  | 1.6 / 0.81mm  | I/A handpiece for 2.7 / 2.8mm incisions            | Curved tip     |  | 5   |  |
| M5561 | 21G     | 0.81mm        | Soft outer I/A handpiece for 2.7 / 2.8mm incisions | 45° angled tip |  | 5   |  |
| M5562 | 21G     | 0.81mm        | Soft outer I/A handpiece for 2.7 / 2.8mm incisions | 60° angled tip |  | 5   |  |
| M5563 | 21G     | 0.81mm        | Soft outer I/A handpiece for 2.7 / 2.8mm incisions | 90° angled tip |  | 5   |  |
| M5571 | 24G     | 0.56mm        | Soft outer I/A handpiece for 2.2mm incisions       | 45° angled tip |  | 5   |  |
| M5572 | 24G     | 0.56mm        | Soft outer I/A handpiece for 2.2mm incisions       | 60° angled tip |  | 5   |  |
| M5573 | 24G     | 0.56mm        | Soft outer I/A handpiece for 2.2mm incisions       | 90° angled tip |  | 5   |  |

# Irrigating/Aspirating Cannulas and Handpieces

For the aspiration of cortex while simultaneously maintaining the Anterior Chamber with continuous irrigating flow.

## Simcoe Irrigating / Aspirating Cannulas

|                      |                                                                                               |                                     |                                                             |                                  |                                                                                     |     |
|----------------------|-----------------------------------------------------------------------------------------------|-------------------------------------|-------------------------------------------------------------|----------------------------------|-------------------------------------------------------------------------------------|-----|
| <b>NEW<br/>M6613</b> | <b>23 G</b>  | <b>23 G / 23 G</b><br>0.64 / 0.64mm | I/A Cannula Standard.<br>200mm silicone tubing.<br>(Simcoe) | 0.4mm port.<br>2.65 OD x 0.65 ID |  | BOX |
|                      |                                                                                               |                                     |                                                             |                                  |                                                                                     | 10  |
| <b>NEW<br/>M6614</b> | <b>23 G</b>  | <b>23 G / 23 G</b><br>0.64 / 0.64mm | I/A Cannula Standard.<br>270mm silicone tubing.<br>(Simcoe) | 0.4mm port.<br>2.65 OD x 0.65 ID |  | BOX |
|                      |                                                                                               |                                     |                                                             |                                  |                                                                                     | 10  |
| <b>NEW<br/>M6617</b> | <b>23 G</b>  | <b>23 G / 23 G</b><br>0.64 / 0.64mm | I/A Cannula Standard.<br>200mm silicone tubing.<br>(Simcoe) | 0.3mm port.<br>2.0 OD x 1.0 ID   |  | BOX |
|                      |                                                                                               |                                     |                                                             |                                  |                                                                                     | 10  |

## Irrigating Microsurgical Handle

|              |                                                                                     |                                                                                       |     |
|--------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|-----|
| <b>M1000</b> |  | Irrigating Microsurgical Handle with female luer taper and a male luer lock connector | BOX |
|              |                                                                                     |                                                                                       | 5   |



# Irrigating/Aspirating Cannulas and Handpieces

For the aspiration of cortex with one handle while simultaneously maintaining the Anterior Chamber with continuous irrigating flow with the other handle.

## Bimanual Cannulas

| Model | Gauge | Image                                                                             | Size                       | Description                         | Tip              | BOX |     |
|-------|-------|-----------------------------------------------------------------------------------|----------------------------|-------------------------------------|------------------|-----|-----|
|       |       |                                                                                   |                            |                                     |                  | 10  | 100 |
| M6524 | 21 G  |  | 21 G x 5/8"<br>0.81 x 16mm | Bimanual aspiration cannula. Curved | 0.3mm top port   |     |     |
| M5405 | 21 G  |  | 21 G x 5/8"<br>0.81 x 16mm | Bimanual aspiration cannula. Curved | 45° bevelled tip | 10  | 100 |

## Coaxial Handpieces

| Model | Gauge | Image                                                                               | Size                       | Description                           | Tip                               | BOX |  |
|-------|-------|-------------------------------------------------------------------------------------|----------------------------|---------------------------------------|-----------------------------------|-----|--|
|       |       |                                                                                     |                            |                                       |                                   | 5   |  |
| M5501 | 23 G  |    | 23 G x 5/8"<br>0.64 x 16mm | Bimanual aspiration handpiece. Curved | 0.4mm side port. Micro-etched tip |     |  |
| M5502 | 23 G  |    | 23 G x 5/8"<br>0.64 x 16mm | Bimanual aspiration handpiece. Curved | 0.3mm side port. Micro-etched tip |     |  |
| M5512 | 21 G  |    | 21 G x 5/8"<br>0.81 x 16mm | Bimanual aspiration handpiece. Curved | 0.3mm side port. Micro-etched tip |     |  |
| M5520 | 21 G  |  | 21 G x 5/8"<br>0.81 x 16mm | Bimanual aspiration handpiece. Curved | 0.4mm side port. Micro-etched tip |     |  |
| M5511 | 20 G  |  | 20 G x 3/4"<br>0.9 x 19mm  | Bimanual aspiration handpiece. Curved | 0.3mm side port. Micro-etched tip |     |  |
| M5518 | 23 G  |  | 23 G x 5/8"<br>0.64 x 16mm | Bimanual aspiration handpiece. Curved | 0.3mm side port.                  |     |  |
| M5515 | 21 G  |  | 21 G x 7/8"<br>0.81 x 22mm | Bimanual aspiration handpiece. Curved | 0.3mm side port.                  |     |  |
| M5519 | 23 G  |  | 23 G x 5/8"<br>0.64 x 16mm | Bimanual irrigation handpiece. Curved | Ultra thin wall. 45° bevelled tip |     |  |
| M5505 | 21 G  |  | 21 G x 5/8"<br>0.81 x 16mm | Bimanual irrigation handpiece. Curved | 45° bevelled tip                  |     |  |
| M5513 | 23 G  |  | 23 G x 5/8"<br>0.64 x 16mm | Bimanual irrigation handpiece. Curved | Thin wall. 2 x 0.4mm side ports   |     |  |
| M5521 | 21G   |  | 21 G x 5/8"<br>0.81 x 16mm | Bimanual irrigation handpiece. Curved | Thin wall. 2 x 0.5mm side ports   |     |  |
| M5514 | 21G   |  | 21 G x 5/8"<br>0.81 x 16mm | Bimanual irrigation handpiece. Curved | 2 x 0.5mm side ports              |     |  |
| M5510 | 20G   |  | 20 G x 3/4"<br>0.9 x 19mm  | Bimanual irrigation handpiece. Curved | 2 x 0.5mm side ports              |     |  |
| M5523 | 23G   |  | 23 G x 5/8"<br>0.64 x 16mm | Bimanual irrigation handpiece. Curved | Open end                          |     |  |


# Irrigating/Aspirating Cannulas and Handpieces

For the aspiration of cortex with one handle while simultaneously maintaining the Anterior Chamber with continuous irrigating flow with the other handle.

## Bimanual I/A Systems

|                    |                          |                                                                                     |                                              |                                                                                       |         |
|--------------------|--------------------------|-------------------------------------------------------------------------------------|----------------------------------------------|---------------------------------------------------------------------------------------|---------|
| <b>M5537 : SET</b> | <b>20 G</b>              |    |                                              |                                                                                       |         |
|                    | <b>21 G</b>              |                                                                                     |                                              |                                                                                       |         |
| <b>Irrigating</b>  | 20 G x ¾"<br>0.9 x 19mm  | Bimanual irrigation handpiece. Curved                                               | 2 x 0.5mm side ports                         |    | BOX     |
| <b>Aspirating</b>  | 21 G x ⅞"<br>0.81 x 22mm | Bimanual aspiration handpiece. Curved                                               | 0.3mm side port                              |    | 5 PAIRS |
| <b>M5531 : SET</b> |                          |                                                                                     |                                              |                                                                                       |         |
| <b>M5531 : SET</b> | <b>21 G</b>              |    |                                              |                                                                                       |         |
|                    |                          |                                                                                     |                                              |                                                                                       |         |
| <b>Irrigating</b>  | 21 G x ⅝"<br>0.81 x 16mm | Bimanual irrigation handpiece. Curved                                               | Thin wall.<br>2 x 0.5mm side ports           |    | BOX     |
| <b>Aspirating</b>  |                          | Bimanual aspiration handpiece. Curved                                               | 0.4mm side port.<br>Micro-etched tip         |    | 5 PAIRS |
| <b>M5331 : SET</b> |                          |                                                                                     |                                              |                                                                                       |         |
| <b>M5331 : SET</b> | <b>21 G</b>              |   |                                              |                                                                                       |         |
|                    |                          |                                                                                     |                                              |                                                                                       |         |
| <b>Irrigating</b>  | 21 G x ⅝"<br>0.81 x 16mm | Bimanual irrigation handpiece. Curved                                               | 2 x 0.4mm side ports.<br>Tapered tip         |  | BOX     |
| <b>Aspirating</b>  |                          | Bimanual aspiration handpiece. Curved                                               | 0.3mm side port.<br>Micro-etched tapered tip |  | 5 PAIRS |
| <b>M5535 : SET</b> |                          |                                                                                     |                                              |                                                                                       |         |
| <b>M5535 : SET</b> | <b>23 G</b>              |  |                                              |                                                                                       |         |
|                    |                          |                                                                                     |                                              |                                                                                       |         |
| <b>Irrigating</b>  | 23 G x ⅝"<br>0.64 x 16mm | Bimanual irrigation handpiece. Curved                                               | Thin wall.<br>2 x 0.4mm side ports           |  | BOX     |
| <b>Aspirating</b>  |                          | Bimanual aspiration handpiece. Curved                                               | 0.3mm side port.<br>Micro-etched tip         |  | 5 PAIRS |
| <b>M5536 : SET</b> |                          |                                                                                     |                                              |                                                                                       |         |
| <b>M5536 : SET</b> | <b>23 G</b>              |  |                                              |                                                                                       |         |
|                    |                          |                                                                                     |                                              |                                                                                       |         |
| <b>Irrigating</b>  | 23 G x ⅝"<br>0.64 x 16mm | Bimanual irrigation handpiece. Curved                                               | Ultra thin wall.<br>45° bevelled tip         |  | BOX     |
| <b>Aspirating</b>  |                          | Bimanual aspiration handpiece. Curved                                               | 0.3mm side port                              |  | 5 PAIRS |
| <b>M5538 : SET</b> |                          |                                                                                     |                                              |                                                                                       |         |
| <b>M5538 : SET</b> | <b>25 G</b>              |  |                                              |                                                                                       |         |
|                    |                          |                                                                                     |                                              |                                                                                       |         |
| <b>Irrigating</b>  | 25 G x ⅝"<br>0.5 x 16mm  | Bimanual irrigation handpiece. Curved                                               | Ultra thin wall.<br>45° bevelled tip         |  | BOX     |
| <b>Aspirating</b>  |                          | Bimanual aspiration handpiece. Curved                                               | 0.3mm side port                              |  | 5 PAIRS |

For more information of other systems available, please contact Sterimedix or your local distributor.



## Lachrymal Cannulas

For probing and irrigating the lachrymal ducts.

|       |                                                                                          |                              |                                   | BOX                             |                                                                                             |
|-------|------------------------------------------------------------------------------------------|------------------------------|-----------------------------------|---------------------------------|---------------------------------------------------------------------------------------------|
| M9275 | 26 G  | 26 G x 3/4"<br>0.45 x 20mm   | Lachrymal cannula. Short straight |                                 |  10100 |
| M9276 | 26 G  | 26 G x 1 1/4"<br>0.45 x 32mm | Lachrymal cannula. Curved         |                                 |  10100 |
| M9277 | 25 G  | 25 G x 1 1/8"<br>0.5 x 29mm  | Lachrymal cannula. Curved         |                                 |  10100 |
| M9274 | 21 G  | 21 G x 1 1/2"<br>0.81 x 38mm | Lachrymal cannula. Curved         |                                 |  10100 |
| M9278 | 25 G  | 25 G x 7/8"<br>0.5 x 22mm    | Lachrymal cannula. Straight       | 0.2mm side port                 |  10100 |
| M9272 | 23 G  | 23 G x 1/2"<br>0.64 x 13mm   | Lachrymal cannula. Straight       | 0.3mm side port.<br>Tapered tip |  10100 |

# Vitreoretinal Products

## Trocar Sets

Provides three entry conduits to the posterior segment of the eye, to facilitate the introduction of various micro surgical ophthalmic instruments .

|              |      |                                                                                     |      |                                                                                                                                        | BOX |
|--------------|------|-------------------------------------------------------------------------------------|------|----------------------------------------------------------------------------------------------------------------------------------------|-----|
|              |      |                                                                                     |      |                                                                                                                                        |     |
| M4510        | 23 G |    | 23 G | One-step trocar set containing:<br>• 3 x self-sealing valved trocars<br>• 1 x 3.5mm scleral fixation plate<br>• 1 x 185mm BSS infusion | 5   |
| M4515        | 25 G |    | 25 G | One-step trocar set containing:<br>• 3 x self-sealing valved trocars<br>• 1 x 3.5mm scleral fixation plate<br>• 1 x 185mm BSS infusion | 5   |
| NEW<br>M4512 | 23 G |    | 23 G | One-step valved trocar cannula<br>& scleral fixation plate                                                                             | 5   |
| NEW<br>M4517 | 25 G |  | 25 G | One-step valved trocar cannula<br>& scleral fixation plate                                                                             | 5   |



# Vitreoretinal Products

For specific infusion and aspiration functions during retinal procedures

## Vitreoretinal Cannulas

| NEW<br>M4486 | 27 G |  | 27 G x 1 1/4"<br>0.4 x 32mm  | Aspiration cannula |                     | BOX |     |
|--------------|------|-----------------------------------------------------------------------------------|------------------------------|--------------------|--------------------------------------------------------------------------------------------------------|-----|-----|
|              |      |                                                                                   |                              |                    |                                                                                                        | 10  | 100 |
| M4477        | 25 G |  | 25 G x 1 1/4"<br>0.5 x 32mm  | Aspiration cannula |                     | 10  | 100 |
| M4478        | 23 G |  | 23 G x 1 1/4"<br>0.64 x 32mm | Aspiration cannula |                     | 10  | 100 |
| M4478B       | 23 G |  | 23 G x 1 1/4"<br>0.64 x 32mm | Aspiration cannula | Ultra thin wall<br> | 10  | 100 |
| M4471        | 20 G |  | 20 G x 1 1/4"<br>0.9 x 32mm  | Aspiration cannula |                     | 10  | 100 |

## Viscous Fluid Infusion Cannulas

|       |      |                                                                                     |                            |                                 |                                                                                                      |     |     |
|-------|------|-------------------------------------------------------------------------------------|----------------------------|---------------------------------|------------------------------------------------------------------------------------------------------|-----|-----|
| M4453 | 25 G |    | 25 G x 1/4"<br>0.5 x 7mm   | Viscous fluid injection cannula |                   | BOX |     |
|       |      |                                                                                     |                            |                                 |                                                                                                      | 10  | 100 |
| M4452 | 23 G |   | 23 G x 1/4"<br>0.64 x 8mm  | Viscous fluid injection cannula | Bevelled tip<br> | 10  | 100 |
| M4447 | 23 G |  | 23 G x 1/25"<br>0.64 x 4mm | Viscous fluid injection cannula |                 | 10  | 100 |
| M4448 | 25 G |  | 25 G x 1/25"<br>0.5 x 4mm  | Viscous fluid injection cannula |                 | 10  | 100 |

## Dual-Bore Cannulas

|        |         |                                                                                     |              |                                                |                                                                                       |     |     |
|--------|---------|-------------------------------------------------------------------------------------|--------------|------------------------------------------------|---------------------------------------------------------------------------------------|-----|-----|
| M4474  | 20/25 G |  | 0.9 / 0.5mm  | Perfluorocarbon coaxial<br>(dual bore) cannula |  | BOX |     |
|        |         |                                                                                     |              |                                                |                                                                                       | 10  | 100 |
| M4499  | 23/30 G |  | 0.64 / 0.3mm | Perfluorocarbon coaxial<br>(dual bore) cannula |  | 10  | 100 |
| M4499A | 25/30 G |  | 0.5 / 0.3mm  | Perfluorocarbon coaxial<br>(dual bore) cannula |  | 10  | 100 |

# Vitreoretinal Products


For specific infusion and aspiration functions during retinal procedures

## Silicone Brush Cannulas

| M4449 | 25 G |  | 25 G x 1 3/8"<br>0.5 x 35mm  | Silicone brush | Silicone extends 0.7mm<br>beyond end of cannula  |  | BOX |     |
|-------|------|-----------------------------------------------------------------------------------|------------------------------|----------------|--------------------------------------------------|-------------------------------------------------------------------------------------|-----|-----|
|       |      |                                                                                   |                              |                |                                                  |                                                                                     | 10  | 100 |
| M4459 | 23 G |  | 23 G x 1 1/4"<br>0.64 x 32mm | Silicone brush | Silicone extends 0.7mm<br>beyond end of cannula  |  | 10  | 100 |
| M4475 | 20 G |  | 20 G x 1 1/4"<br>0.9 x 32mm  | Silicone brush | Silicone extends 1.25mm<br>beyond end of cannula |  | 10  | 100 |

## Soft Tip Cannulas

| M4456 | 25 G |  | 25 G x 1 1/4"<br>0.5 x 32mm  | Soft tip | Silicone extends 1mm<br>beyond end of cannula |  | BOX |     |
|-------|------|-----------------------------------------------------------------------------------|------------------------------|----------|-----------------------------------------------|-------------------------------------------------------------------------------------|-----|-----|
|       |      |                                                                                   |                              |          |                                               |                                                                                     | 10  | 100 |
| M4458 | 23 G |  | 23 G x 1 1/4"<br>0.64 x 32mm | Soft tip | Silicone extends 2mm<br>beyond end of cannula |  | 10  | 100 |
| M4457 | 20 G |  | 20 G x 1 1/4"<br>0.9 x 32mm  | Soft tip | Silicone extends 1mm<br>beyond end of cannula |  | 10  | 100 |




# Vitreoretinal Products

For specific infusion and aspiration functions during retinal procedures

## Backflush Handles

|              |                                                                                   |                                                                                               |     |
|--------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|-----|
| <b>M4489</b> |  | Backflush handle for passive aspiration with finger control and a luer lok cannula connection | BOX |
|              |                                                                                   |                                                                                               | 5   |


## 'Active / Passive' Backflush Handles and Combination Sets

|                        |                                                                                                 |                              |                                                                   |                                                                                       |          |
|------------------------|-------------------------------------------------------------------------------------------------|------------------------------|-------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------|
| <b>M4503 : SET</b>     | <b>25 G</b>    | 25 G x 1 1/4"<br>0.5 x 32mm  | 'Active/passive' backflush handle with 25G soft tip cannula       |    | BOX<br>5 |
| <b>M4502 : SET</b>     | <b>23 G</b>    | 23 G x 1 1/4"<br>0.64 x 32mm | 'Active/passive' backflush handle with 23G soft tip cannula       |    | 5        |
| <b>M4501 : SET</b>     | <b>20 G</b>    | 20 G x 1 1/4"<br>0.9 x 32mm  | 'Active/passive' backflush handle with 20G soft tip cannula       |    | 5        |
| <b>NEW M4505 : SET</b> | <b>25 G</b>    | 25 G x 1 1/4"<br>0.5 x 32mm  | 'Active/passive' backflush handle with 25G silicone brush cannula |    | 5        |
| <b>M4492 : SET</b>     | <b>23 G</b>    | 23 G x 1 1/4"<br>0.64 x 32mm | 'Active/passive' backflush handle with 23G silicone brush cannula |    | 5        |
| <b>M4500 : SET</b>     | <b>25 G</b>  | 25 G x 1 1/4"<br>0.5 x 32mm  | 'Active/passive' backflush handle with 25G aspiration cannula     |   | 5        |
| <b>M4497 : SET</b>     | <b>23 G</b>  | 23 G x 1 1/4"<br>0.64 x 32mm | 'Active/passive' backflush handle with 23G aspiration cannula     |  | 5        |
| <b>M4496 : SET</b>     | <b>20 G</b>  | 20 G x 1 1/4"<br>0.9 x 32mm  | 'Active/passive' backflush handle with 20G aspiration cannula     |  | 5        |

## Backflush Handles

|               |                                                                                     |                                                                                                             |     |
|---------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|-----|
| <b>M4489A</b> |  | 'Active/passive' backflush handle with finger control. Allows conversion from active to passive aspiration. | BOX |
|               |                                                                                     |                                                                                                             | 5   |

## 'Active / Passive' Backflush Handles and Combination Sets

|                        |                                                                                                 |                              |                                                                   |                                                                                       |          |
|------------------------|-------------------------------------------------------------------------------------------------|------------------------------|-------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------|
| <b>M4504 : SET</b>     | <b>25 G</b>  | 25 G x 1 1/4"<br>0.5 x 32mm  | 'Active/passive' backflush handle with 25G soft tip cannula       |  | BOX<br>5 |
| <b>M4530 : SET</b>     | <b>23 G</b>  | 23 G x 1 1/4"<br>0.64 x 32mm | 'Active/passive' backflush handle with 23G soft tip cannula       |  | 5        |
| <b>M4531 : SET</b>     | <b>20 G</b>  | 20 G x 1 1/4"<br>0.9 x 32mm  | 'Active/passive' backflush handle with 20G soft tip cannula       |  | 5        |
| <b>NEW M4507 : SET</b> | <b>25 G</b>  | 25 G x 1 1/4"<br>0.5 x 32mm  | 'Active/passive' backflush handle with 25G silicone brush cannula |  | 5        |
| <b>NEW M4509 : SET</b> | <b>23 G</b>  | 23 G x 1 1/4"<br>0.64 x 32mm | 'Active/passive' backflush handle with 23G silicone brush cannula |  | 5        |

# Vitreoretinal Products

For specific infusion and aspiration functions during retinal procedures

## Vitreoretinal Heavy Liquid Infusion Handles

|       |                                                                                   |                                                                            |     |
|-------|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------|-----|
| M4498 |  | Heavy liquid infusion handle with finger control and a luer lok connection | BOX |
|       |                                                                                   |                                                                            | 5   |

## Charles Flute Handles

|       |                                                                                   |                                                                                                             |     |
|-------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|-----|
| M4493 |  | 'Active/passive' backflush handle with finger control. Allows conversion from active to passive aspiration. | BOX |
|       |                                                                                   |                                                                                                             | 5   |

## Charles Flute Handles and Combination Sets

|             |                                                                                        |                              |                                                    |                                                                                     |     |
|-------------|----------------------------------------------------------------------------------------|------------------------------|----------------------------------------------------|-------------------------------------------------------------------------------------|-----|
| M4495 : SET | 23 G  | 23 G x 1 1/4"<br>0.64 x 32mm | 'Charles' flute handle with 23G aspiration cannula |  | BOX |
|             |                                                                                        |                              |                                                    |                                                                                     | 5   |



Sterimedix products are packed sterile in a single blister in space saving boxes designed to minimise waste.

## Protector for Kit Packers

A new labelled protector for bulk non-sterile products, to provide increased flexibility and reduce packaging waste.

Designed to protect the cannulas and needle tips by ensuring that no distortion or deformation can occur to formed or angled products such as cystotomes.

For further information, please contact us at: [sales@sterimedix](mailto:sales@sterimedix).



[www.sterimedix.com](http://www.sterimedix.com)

| PRODUCT CODE | SIZE          | DESCRIPTION                                    | PAGE | NO. PER BOX |     |
|--------------|---------------|------------------------------------------------|------|-------------|-----|
| M0275        | 23 G x 1 1/2" | Retrobulbar. Tri-facet                         | 5    | 10          | 100 |
| M0278        | 19 G x 1"     | Sub-tenon cannula                              | 5    | 10          | 100 |
| M0279        | 20 G x 1"     | Sub-tenon cannula                              | 5    | 10          | 100 |
| M0637        | 25 G x 1 1/2" | Retrobulbar (Atkinson)                         | 5    | 10          | 100 |
| M0638        | 23 G x 1 1/2" | Retrobulbar (Atkinson)                         | 5    | 10          | 100 |
| M0641        | 27 G x 7/8"   | Retrobulbar (Atkinson)                         | 5    | 10          | 100 |
| M0641A       | 27 G x 1 1/4" | Retrobulbar (Atkinson)                         | 5    | 10          | 100 |
| M0642        | 25 G x 1 1/4" | Retrobulbar (Atkinson)                         | 5    | 10          | 100 |
| M0642A       | 25 G x 7/8"   | Retrobulbar (Atkinson)                         | 5    | 10          | 100 |
| M0643        | 23 G x 1 1/4" | Retrobulbar (Atkinson)                         | 5    | 10          | 100 |
| M1000        | 100mm         | Irrigating microsurgical handle                | 15   | 5           |     |
| M2270        | 30 G x 7/8"   | Anterior chamber cannula. Angled 35° (Rycroft) | 6    | 10          | 100 |
| M2270A       | 30 G x 7/8"   | Anterior chamber cannula. Angled 35° (Rycroft) | 6    | 10          | 100 |
| M2272        | 26 G x 7/8"   | Viscoelastic cannula. Angled 35°               | 7    | 10          | 100 |
| M2273        | 30 G x 7/8"   | Anterior chamber cannula. Angled 35° (Rycroft) | 6    | 10          | 100 |
| M2273A       | 27 G x 7/8"   | Anterior chamber cannula. Angled 35° (Rycroft) | 6    | 10          | 100 |
| M2273B       | 20 G x 7/8"   | Anterior chamber cannula. Angled 45° (Rycroft) | 6    | 10          | 100 |
| M2273C       | 27 G x 7/8"   | Hydrodissector. Curved (Sauter)                | 10   | 10          | 100 |
| M2273D       | 25 G x 7/8"   | Hydrodissector. Curved (Sauter)                | 10   | 10          | 100 |
| M2273E       | 27 G x 7/8"   | Hydrodissector. Curved (Helsinki)              | 10   | 10          | 100 |
| M2273F       | 25 G x 7/8"   | Hydrodissector. Curved (Helsinki)              | 10   | 10          | 100 |
| M2273G       | 23 G x 7/8"   | Anterior chamber cannula. Angled 35° (Rycroft) | 6    | 10          | 100 |
| M2273H       | 27 G x 7/8"   | Hydrodissector. Curved                         | 10   | 10          | 100 |
| M2274        | 25 G x 7/8"   | Hydrodissector. Curved (Kellan)                | 10   | 10          | 100 |
| M2275        | 25 G x 7/8"   | Anterior chamber cannula. Angled 35° (Rycroft) | 6    | 10          | 100 |
| M2279        | 27 G x 7/8"   | Anterior chamber cannula. Angled 45°           | 6    | 10          | 100 |
| M2280        | 27 G x 7/8"   | Transzonular injection cannula. Angled 33°     | 6    | 10          | 100 |
| M2639        | 25 G x 1"     | Anterior chamber cannula. Angled 35° (Rycroft) | 6    | 10          | 100 |
| M2639A       | 25 G x 7/8"   | Viscoelastic cannula. Angled 35°               | 7    | 10          | 100 |
| M2642        | 21 G x 7/8"   | Anterior chamber cannula. Angled 35° (Rycroft) | 6    | 10          | 100 |
| M2657        | 20 G x 7/8"   | Anterior chamber maintainer (Lewicky)          | 6    | 10          |     |
| M2658        | 23 G x 9/64"  | Anterior chamber maintainer (Lewicky)          | 6    | 10          |     |
| M2801        | 25 G x 7/8"   | Irrigating lasik cannula. 2 ports (Lin)        | 4    | 10          | 100 |
| M2804        | 25 G x 7/8"   | Irrigating lasik cannula. 6 ports              | 4    | 10          | 100 |
| M2817        | 27 G x 7/8"   | Irrigating lasik cannula. 4 ports              | 4    | 10          | 100 |
| M2822        | 25 G x 7/8"   | Irrigating lasik cannula. 6 ports              | 4    | 10          | 100 |
| M2824        | 26 G x 7/8"   | Irrigating lasik cannula (Slade)               | 4    | 10          | 100 |
| M3000D       | 25 G x 5/8"   | Cystotome. Short curve. Formed                 | 9    | 10          | 100 |
| M3000E       | 27 G x 5/8"   | Cystotome. Short curve. Formed                 | 9    | 10          | 100 |
| M3000F       | 27 G x 5/8"   | Cystotome. Formed                              | 9    | 10          | 100 |
| M3000G       | 25 G x 5/8"   | Cystotome. Formed                              | 9    | 10          | 100 |
| M3001 : NEW  | 27 G x 5/8"   | Endocapsular cystotome. Formed (Joyce)         | 9    | 10          | 100 |
| M3602        | 25 G x 5/8"   | Cystotome. Side cutting. Formed                | 8    | 10          | 100 |
| M3610        | 25 G x 5/8"   | Cystotome. Formed                              | 8    | 10          | 100 |
| M3610B       | 27 G x 1 1/2" | Cystotome. Formed                              | 8    | 10          | 100 |
| M3610C       | 25 G x 5/8"   | Cystotome. Formed. Angled 35°                  | 8    | 10          | 100 |
| M3610S       | 25 G x 5/8"   | Cystotome. Short tip. Formed                   | 8    | 10          | 100 |
| M3615        | 27 G x 7/8"   | Long cystotome. Formed                         | 8    | 10          | 100 |
| M3626        | 27 G x 5/8"   | Cystotome. Angled                              | 8    | 10          | 100 |
| M3626A       | 25 G x 5/8"   | Cystotome. Angled                              | 8    | 10          | 100 |
| M3627        | 25 G x 5/8"   | Cystotome. Short curve. Formed                 | 8    | 10          | 100 |
| M3627A       | 27 G x 5/8"   | Cystotome. Short curve. Formed                 | 8    | 10          | 100 |

| PRODUCT CODE | SIZE                  | DESCRIPTION                                            | PAGE | NO. PER BOX |     |
|--------------|-----------------------|--------------------------------------------------------|------|-------------|-----|
| M3629        | 25 G x 5/8"           | Cystotome. Reverse tip. Formed                         | 8    | 10          | 100 |
| M3630        | 30 G x 5/8"           | Cystotome. Formed                                      | 8    | 10          | 100 |
| M3899        | 25 G x 7/8"           | Viscoexpression cannula. Flat tip. Angled (Corydon)    | 11   | 10          | 100 |
| M3899A       | 25 G x 7/8"           | Viscoexpression cannula. Rounded tip. Angled (Corydon) | 11   | 10          | 100 |
| M3900        | 25 G x 7/8"           | Hydrodissector. Angled. Tapered and flattened tip      | 10   | 10          | 100 |
| M3902        | 27 G x 7/8"           | Hydrodissector. Angled. Tapered and flattened tip      | 10   | 10          | 100 |
| M3904        | 27 G x 3/4"           | Hydrodissector. Angled 90° (Chang)                     | 11   | 10          | 100 |
| M3906A       | 25 G x 7/8"           | Hydrodissector. Rounded 'J'-shaped tip                 | 11   | 10          | 100 |
| M3907        | 25 G x 7/8"           | Hydrodissector. Angled 45°. Tapered to 0.3mm           | 11   | 10          | 100 |
| M3908        | 25 G x 7/8"           | Hydrodissector. Angled 45°. Tapered to 0.3mm           | 11   | 10          | 100 |
| M3909        | 23 G x 7/8"           | Viscoelastic cannula. Angled 45°                       | 7    | 10          | 100 |
| M3910        | 21 G x 7/8"           | Viscoelastic cannula. Straight                         | 7    | 10          | 100 |
| M3910A       | 21 G x 7/8"           | Viscoelastic cannula. Angled 45°                       | 7    | 10          | 100 |
| M3911        | 27 G x 7/8"           | Hydrodissector. Angled. Flattened tip                  | 11   | 10          | 100 |
| M3911A       | 27 G x 7/8"           | Hydrodissector. Angled. Round tip                      | 11   | 10          | 100 |
| M3912        | 23 G x 7/8"           | Hydrodissector. Angled 30° (Mackool)                   | 11   | 10          | 100 |
| M3914        | 23 G x 7/8"           | Irrigating lasik cannula. Formed. Flattened tip        | 4    | 10          | 100 |
| M3917        | 27 G x 7/8"           | Viscoelastic cannula. Angled. Formed                   | 7    | 10          | 100 |
| M3919        | 27 G x 7/8"           | Hydrodissector. Tapered tip (Akahoshi)                 | 10   | 10          | 100 |
| M3920 : NEW  | 27 G x 7/8"           | Hydrodissector. Curved. (Frankfurt)                    | 11   | 10          | 100 |
| M3921 : NEW  | 25 G x 7/8"           | Hydrodissector. Curved and tapered                     | 11   | 10          | 100 |
| M3922 : NEW  | 27 G x 7/8"           | Hydrodissector. Formed                                 | 11   | 10          | 100 |
| M3926A       | 27 G x 7/8"           | Viscoelastic cannula. Angled 35°                       | 7    | 10          | 100 |
| M3927        | 27 G x 7/8"           | Viscoelastic cannula. Angled 25°                       | 7    | 10          | 100 |
| M3992        | 30 G x 7/8"           | Hydrodissector. Angled. Flattened tip                  | 10   | 10          | 100 |
| M4447        | 23 G x 1/4"           | Viscous fluid injection cannula                        | 20   | 10          | 100 |
| M4448        | 25 G x 1/4"           | Viscous fluid injection cannula                        | 20   | 10          | 100 |
| M4449        | 25 G x 1 3/8"         | Silicone brush. 0.7mm                                  | 21   | 10          | 100 |
| M4452        | 23 G x 1/4"           | Viscous fluid injection cannula. Bevelled tip          | 20   | 10          | 100 |
| M4453        | 25 G x 1/4"           | Viscous fluid injection cannula                        | 20   | 10          | 100 |
| M4456        | 25 G x 1 1/4"         | Soft tip cannula. 1mm                                  | 21   | 10          | 100 |
| M4457        | 20 G x 1 1/4"         | Soft tip cannula. 1mm                                  | 21   | 10          | 100 |
| M4458        | 23 G x 1 1/4"         | Soft tip cannula. 2mm                                  | 21   | 10          | 100 |
| M4459        | 23 G x 1 1/4"         | Silicone brush. 0.7mm                                  | 21   | 10          | 100 |
| M4471        | 20 G x 1 1/4"         | Aspiration cannula                                     | 20   | 10          | 100 |
| M4474        | 20 G / 25 G           | Perfluorocarbon coaxial cannula (dual bore)            | 20   | 10          | 100 |
| M4475        | 20 G x 1 1/4"         | Silicone brush. 1.25mm                                 | 21   | 10          | 100 |
| M4477        | 25 G x 1 1/4"         | Aspiration cannula                                     | 20   | 10          | 100 |
| M4478        | 23 G x 1 1/4"         | Aspiration cannula                                     | 20   | 10          | 100 |
| M4478B       | 23 G x 1 1/4"         | Aspiration cannula. Ultra thin wall                    | 20   | 10          | 100 |
| M4486 : NEW  | 27 G x 1 1/4"         | Aspiration cannula                                     | 20   | 10          | 100 |
| M4489        | 100mm                 | Backflush handle                                       | 22   | 5           |     |
| M4489A       | 100mm                 | Active / passive backflush handle                      | 22   | 5           |     |
| M4492 : SET  | 100mm / 23 G x 1 1/4" | Backflush handle + 23 G silicone brush                 | 22   | 5           |     |
| M4493        | 110mm                 | 'Charles' flute handle                                 | 23   | 5           |     |
| M4495 : SET  | 110mm / 23 G x 1 1/4" | 'Charles' flute handle + 23 G aspiration cannula       | 23   | 5           |     |
| M4496 : SET  | 100mm / 20 G x 1 1/4" | Backflush handle + 20 G aspiration cannula             | 22   | 5           |     |
| M4497 : SET  | 100mm / 23 G x 1 1/4" | Backflush handle + 23 G aspiration cannula             | 22   | 5           |     |
| M4498        | 110mm                 | Heavy liquid infusion handle                           | 23   | 5           |     |
| M4499        | 23 G / 30 G           | Perfluorocarbon coaxial cannula (dual bore)            | 20   | 10          | 100 |
| M4499A       | 25 G / 30 G           | Perfluorocarbon coaxial cannula (dual bore)            | 20   | 10          | 100 |


| PRODUCT CODE | SIZE               | DESCRIPTION                                                                      | PAGE | NO. PER BOX |     |
|--------------|--------------------|----------------------------------------------------------------------------------|------|-------------|-----|
| M4500 : SET  | 100mm / 25 G x 1¼" | Backflush handle + 25 G aspiration cannula                                       | 22   | 5           |     |
| M4501 : SET  | 100mm / 20 G x 1¼" | Backflush handle + 20 G soft tip cannula                                         | 22   | 5           |     |
| M4502 : SET  | 100mm / 23 G x 1¼" | Backflush handle + 23 G soft tip cannula                                         | 22   | 5           |     |
| M4503 : SET  | 100mm / 25 G x 1¼" | Backflush handle + 25 G soft tip cannula                                         | 22   | 5           |     |
| M4504 : SET  | 100mm / 25 G x 1¼" | Active / passive backflush handle + 25 G soft tip cannula                        | 22   | 5           |     |
| M4505 : SET  | 100mm / 25 G x 1¼" | Backflush handle + 25 G silicone brush                                           | 22   | 5           |     |
| M4507 : SET  | 100mm / 25 G x 1¼" | Active / passive backflush handle + 25 G silicone brush                          | 22   | 5           |     |
| M4509 : SET  | 100mm / 23 G x 1¼" | Active / passive backflush handle + 23 G silicone brush                          | 22   | 5           |     |
| M4510 : SET  | 23 G               | One-step trocar system                                                           | 19   | 5           |     |
| M4512 : SET  | 23 G               | One-step valved trocar cannula + scleral fixation plate                          | 19   | 5           |     |
| M4515 : SET  | 25 G               | One-step trocar system                                                           | 19   | 5           |     |
| M4517 : SET  | 25 G               | One-step valved trocar cannula + scleral fixation plate                          | 19   | 5           |     |
| M4530 : SET  | 100mm / 23 G x 1¼" | Active / passive backflush handle + 23 G silicone brush                          | 22   | 5           |     |
| M4531 : SET  | 100mm / 20 G x 1¼" | Active / passive backflush handle + 20 G silicone brush                          | 22   | 5           |     |
| M4619        | 25 G x 1½"         | Irrigating vectis. 5mm. Single port                                              | 12   | 10          | 100 |
| M4620        | 25 G x 1½"         | Irrigating vectis. 5mm. 2 ports                                                  | 12   | 10          | 100 |
| M5331 : SET  | 21 G x ⅝"          | Bimanual irrigation handpiece. Curved. 2 x 0.4mm side ports. Tapered tip         | 17   | 5 PAIRS     |     |
|              | 21 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.3mm side port. Micro-etched tapered tip |      |             |     |
| M5405        | 21 G x ⅝"          | Bimanual irrigation cannula. Curved. 45° bevelled tip                            | 16   | 10          | 100 |
| M5501        | 23 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.4mm side port. Micro-etched tip         | 16   | 5           |     |
| M5502        | 23 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.3mm side port. Micro-etched tip         | 16   | 5           |     |
| M5505        | 21 G x ⅝"          | Bimanual irrigation handpiece. Curved. Bevelled tip                              | 16   | 5           |     |
| M5510        | 20 G x ¾"          | Bimanual irrigation handpiece. Curved. 2 x 0.5mm side ports                      | 16   | 5           |     |
| M5511        | 20 G x ¾"          | Bimanual aspiration handpiece. Curved. 0.3mm side port. Micro-etched tip         | 16   | 5           |     |
| M5512        | 21 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.3mm side port. Micro-etched tip         | 16   | 5           |     |
| M5513        | 23 G x ⅝"          | Bimanual irrigation handpiece. Thin wall. Curved. 2 x 0.4mm side ports           | 16   | 5           |     |
| M5514        | 21 G x ⅝"          | Bimanual irrigation handpiece. Curved. 2 x 0.5mm side ports                      | 16   | 5           |     |
| M5515        | 21 G x ⅞"          | Bimanual aspiration handpiece. Curved. 0.3mm side port                           | 16   | 5           |     |
| M5518        | 23 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.3mm side port                           | 16   | 5           |     |
| M5519        | 23 G x ⅝"          | Bimanual irrigation handpiece. Ultra thin wall. Curved. Bevelled tip             | 16   | 5           |     |
| M5520        | 21 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.4mm side port. Micro-etched tip         | 16   | 5           |     |
| M5521        | 21 G x ⅝"          | Bimanual irrigation handpiece. Thin wall. Curved. 2 x 0.5mm side ports           | 16   | 5           |     |
| M5523        | 23 G x ⅝"          | Bimanual irrigation handpiece. Ultra thin wall. Curved. Open end                 | 16   | 5           |     |
| M5531        | 21 G x ⅝"          | Bimanual irrigation handpiece. Thin wall. Curved. 2 x 0.5mm side ports           | 17   | 5 PAIRS     |     |
|              | 21 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.4mm side port. Micro-etched tip         |      |             |     |
| M5535        | 23 G x ⅝"          | Bimanual irrigation handpiece. Thin wall. Curved. 2 x 0.4mm side ports           | 17   | 5 PAIRS     |     |
|              | 23 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.3mm side port. Micro-etched tip         |      |             |     |
| M5536        | 23 G x ⅝"          | Bimanual irrigation handpiece. Ultra thin wall. Curved. 45° bevelled tip         | 17   | 5 PAIRS     |     |
|              | 23 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.3mm side port                           |      |             |     |
| M5537        | 20 G x ¾"          | Bimanual irrigation handpiece. Curved. 2 x 0.5mm side ports                      | 17   | 5 PAIRS     |     |
|              | 21 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.3mm side port                           |      |             |     |
| M5538        | 25 G x ⅝"          | Bimanual irrigation handpiece. Ultra thin wall. Curved. 45° bevelled tip         | 17   | 5 PAIRS     |     |
|              | 25 G x ⅝"          | Bimanual aspiration handpiece. Curved. 0.3mm side port                           |      |             |     |
| M5561        | 21 G               | Soft outer I/A handpiece. 45° angled tip                                         | 14   | 5           |     |
| M5562        | 21 G               | Soft outer I/A handpiece. 60° angled tip                                         | 14   | 5           |     |
| M5563        | 21 G               | Soft outer I/A handpiece. 90° angled tip                                         | 14   | 5           |     |
| M5571        | 24 G               | Soft outer I/A handpiece for 2.2mm incision. 45° angled tip                      | 14   | 5           |     |
| M5572        | 24 G               | Soft outer I/A handpiece for 2.2mm incision. 60° angled tip                      | 14   | 5           |     |
| M5573        | 24 G               | Soft outer I/A handpiece for 2.2mm incision. 90° angled tip                      | 14   | 5           |     |
| M5580        | 16 G / 21 G        | I/A handpiece. Straight tip                                                      | 14   | 5           |     |
| M5581        | 16 G / 21 G        | I/A handpiece. Curved tip                                                        | 14   | 5           |     |
| M5582        | 16 G / 21 G        | I/A handpiece. 45° angled tip                                                    | 14   | 5           |     |

| PRODUCT CODE | SIZE          | DESCRIPTION                                                                        | PAGE | NO. PER BOX |     |
|--------------|---------------|------------------------------------------------------------------------------------|------|-------------|-----|
| M5583        | 16 G / 21 G   | I/A handpiece. 60° angled tip                                                      | 14   | 5           |     |
| M5584        | 16 G / 21 G   | I/A handpiece. 90° angled tip                                                      | 14   | 5           |     |
| M5591        | 17 G / 24 G   | I/A handpiece for 2.2mm incision. Curved tip                                       | 14   | 5           |     |
| M5592        | 17 G / 24 G   | I/A handpiece for 2.2mm incision. 45° angled tip                                   | 14   | 5           |     |
| M5593        | 17 G / 24 G   | I/A handpiece for 2.2mm incision. 60° angled tip                                   | 14   | 5           |     |
| M5594        | 17 G / 24 G   | I/A handpiece for 2.2mm incision. 90° angled tip                                   | 14   | 5           |     |
| M6524        | 21 G x 5/8"   | Bimanual aspirating cannula. Curved. 0.3mm top port                                | 16   | 10          | 100 |
| M6575        | 23 G x 5/8"   | Cortex extractor. 0.3mm side port                                                  | 12   | 10          | 100 |
| M6575A       | 23 G x 5/8"   | Cortex extractor. 0.4mm side port                                                  | 12   | 10          | 100 |
| M6575B : NEW | 23 G x 5/8"   | Cortex extractor. Curved. 0.3mm side port                                          | 12   | 10          | 100 |
| M6576        | 23 G x 7/8"   | Cortex extractor. 0.3mm side port                                                  | 12   | 10          | 100 |
| M6576A       | 23 G x 1/2"   | Cortex extractor. 0.3mm side port                                                  | 12   | 10          | 100 |
| M6594        | 23 G x 5/8"   | Cortex extractor. Bevelled end port (Charleux)                                     | 12   | 10          | 100 |
| M6595        | 23 G x 1/2"   | Cortex extractor. Bevelled end port (Charleux)                                     | 12   | 10          | 100 |
| M6596        | 24 G x 1/2"   | Cortex extractor. Bevelled end port (Charleux)                                     | 12   | 10          | 100 |
| M6611        | 21 G x 1 1/2" | I/A cannula. Straight                                                              | 14   | 10          | 100 |
| M6613 : NEW  | 23 G / 23 G   | Standard I/A cannula. 0.4mm port. 200mm silicone tubing 2.65 OD x 0.65 ID (Simcoe) | 15   | 10          |     |
| M6614 : NEW  | 23 G / 23 G   | Standard I/A cannula. 0.4mm port. 270mm silicone tubing 2.65 OD x 0.65 ID (Simcoe) | 15   | 10          |     |
| M6617 : NEW  | 23 G / 23 G   | Standard I/A cannula. 0.3mm port. 200mm silicone tubing 2 OD x 1 ID (Simcoe)       | 15   | 10          |     |
| M7500        | 3.5mm / 4.0mm | Scleral marker. 80mm                                                               | 4    | 10          |     |
| M7510        | 3.5mm         | Scleral fixation / incision template. 75mm                                         | 4    | 10          |     |
| M8604        | 21 G x 1"     | Capsule polisher. Angled. 0.3mm side port. Micro-etched tip                        | 13   | 10          | 100 |
| M8604A       | 23 G x 7/8"   | Capsule polisher. Angled. 0.3mm side port. Micro-etched tip                        | 13   | 10          | 100 |
| M8604B       | 21 G x 7/8"   | Capsule polisher. Angled. 0.3mm side port. Micro-etched tip                        | 13   | 10          | 100 |
| M8605        | 27 G x 7/8"   | Silicone tipped polisher. Angled                                                   | 13   | 10          | 100 |
| M8605A       | 23 G x 7/8"   | Silicone tipped polisher. Angled                                                   | 13   | 10          | 100 |
| M8620C       | 23 G x 1"     | Olive tipped polisher. Curved                                                      | 13   | 10          | 100 |
| M8621        | 23 G x 1"     | Olive tipped polisher. Curved. Micro-etched tip                                    | 13   | 10          | 100 |
| M8622        | 25 G x 1"     | Olive tipped polisher. Angled. Micro-etched tip                                    | 13   | 10          | 100 |
| M8645        | 23 G x 7/8"   | Capsule scraper. Angled. Micro-etched distal tip                                   | 13   | 10          | 100 |
| M9272        | 23 G x 1/2"   | Lachrymal cannula. Straight. 0.3mm side port. Tapered tip                          | 18   | 10          | 100 |
| M9274        | 23 G x 1 1/2" | Lachrymal cannula. Curved                                                          | 18   | 10          | 100 |
| M9275        | 26 G x 3/4"   | Lachrymal cannula. Short straight                                                  | 18   | 10          | 100 |
| M9276        | 26 G x 1 1/4" | Lachrymal cannula. Curved                                                          | 18   | 10          | 100 |
| M9277        | 25 G x 1 1/8" | Lachrymal cannula. Curved                                                          | 18   | 10          | 100 |
| M9278        | 25 G x 7/8"   | Lachrymal cannula. Straight. 0.2mm side port                                       | 18   | 10          | 100 |
| M9280        | 19 G x 1"     | Washout cannula. Angled (Bishop Hamon)                                             | 6    | 10          | 100 |
| M9281        | 16 G x 1"     | Washout cannula. Angled. Flattened tip                                             | 6    | 10          | 100 |

**Sterimedix**<sup>®</sup>  
ophthalmic

# EXPERTISE

in single-use ophthalmic surgical products



Distributed by:

**Sterimedix Limited**

Thornhill Road, North Moons Moat, Redditch, Worcestershire B98 9ND, UK  
T : +44 (0)1527 501480 F : +44 (0)1527 501491 e : sales@sterimedix.com

[www.sterimedix.com](http://www.sterimedix.com)

**Sterimedix**<sup>®</sup> is a registered trademark of Sterimedix Limited © Sterimedix Limited 2019. All rights reserved

Sterimedix policy is one of continuous improvement and we therefore reserve the right to alter specifications and prices without prior notice. All products are sold subject to the company's conditions of sale.