



NYLON

NYLON

Non-absorbable, synthetic surgical monofilament thread




- Polyamide
- Uncoated
- Thread may be blue or black or undyed
- Size 10/0 – 5 USP
- Features
 - Uniform structure with smooth surface
 - Flexible and elastic tightening of knots
 - Very good bindings properties
 - High durability of knot
 - Tensile strength
- Indications
 - Soft tissues (intra-dermal and dermal suturing)
 - General, Pediatric, Plastic, Vascular surgery
 - Urology
 - Ophthalmology
- Period of validity 5 years

yavo

YAVO Sp. z o. o. 97-400 Bełchatów, ul. Bawełniana 17, Polska
tel.: +48 44 632 94 11, fax:+48 44 632 49 17,
e-mail: info@yavo.com.pl, www.yavo.com.pl



NYLON



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



spatula (lanced needle)

length	needles			thread length (cm) color	10/0	9/0	8/0	7/0	6/0	5/0	4/0
	curvature	diameter									
6 mm	3/8	150 µm		7 black	LZ63PA						
				15 black	LZ55PA						
				20 black	LZ37PA						
6 mm micropoint	3/8			30 black	LZ60PA						
6 mm micropoint	3/8	150 µm		30 black	LZ52PA	LZ54PA					
6 mm	3/8	150 µm		30 black	LZ10PA	LZ11PA					
6 mm	3/8	200 µm		30 black	LZ12PA	LZ13PA	LZ26PA	LZ76PA			
				45 black			LZ57PA				
6 mm	3/8	200 µm		30 black	LZ14PA	LZ15PA					
				45 black	LZ61PA						
6,19 mm	3/8	150 µm		30 black	LZ49PA	LZ56PA					
6,19 mm	3/8	200 µm		30 black	LZ18PA	LZ19PA					
6,19 mm	3/8	200 µm		30 black	LZ20PA	LZ21PA					
6,2 mm	3/8			30 black	LZ30PA	LZ31PA					
6,2 mm micropoint	3/8	200 µm		30 black	LZ78PA						
6,4 mm	3/8			30 black			LZ41PA				
6,4 mm	3/8	200 µm		30 black	LZ22PA	LZ23PA	LZ24PA				
6,5 mm	3/8	150 µm		30 black		LZ67PA					
6,5 mm	3/8	220 µm		30 black	LZ35PA	LZ36PA					
6,5 mm micropoint	3/8	220 µm		30 black	LZ64PA	LZ65PA					
				45 black	LZ79PA						

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



spatula (lanced needle)

length	needles		curvature	diameter	thread length (cm)	color	10/0	9/0	8/0	7/0	6/0	5/0	4/0
	curvature	diameter											
6,5 mm micropoint	3/8	150 µm			30	black	LZ69PA						
					45	black			LZ77PA				
8 mm	3/8				30	black					LZ40PA		
5,5 mm	1/2				30	black	LZ28PA						
5,5 mm micropoint	1/2	150 µm			30	black	LZ53PA						
5,6 mm	1/2	150 µm			20	black	LZ03PA	LZ04PA					
5,6 mm	1/2	150 µm			30	black	LZ06PA	LZ07PA	LZ08PA				
6 mm	1/2	150 µm			15	black	LZ33PA						
					30	black	LZ16PA		LZ59PA				
6 mm	1/2	150 µm			30	black	LZ17PA		LZ39PA				
6 mm micropoint	1/2				30	black			LZ70PA				
6,5 mm micropoint	1/2	220 µm			30	black	LZ74PA						
8 mm	1/4	280 µm			45	blue						LZ75PA	

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle

length	needles curvature		thread length (cm) color	10/0	9/0	8/0	7/0	6/0	5/0	4/0
4 mm	3/8		13 black	RS115PA	RS248PA					
			15 black	RS77PA						
5 mm	3/8		13 black	RS212PA	RS214PA	RS224PA				
			15 black	RS130PA						
			45 black	RS200PA						
6,4 mm	3/8		13 black			RS245PA				
			15 black		RS02PA	RS03PA				
6,4 mm	3/8		45 black	RS89PA	RS04PA	RS05PA				
8 mm	3/8		60 blue				RS182PA			
9,5 mm	3/8		75 blue				RS08PA	RS09PA		
10 mm	3/8		45 blue			RS201PA		RS10PA		
11 mm	3/8		60 blue				RS106PA			
12 mm	3/8		45 blue					RS11PA	RS12PA	
12 mm	3/8		45 blue					RS79PA		
13 mm	3/8		60 blue					RS13PA		
			75 blue			RS158PA	RS159PA	RS184PA	RS185PA	
18 mm	3/8		90 blue						RS160PA	

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle

length	needles curvature	thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
			19 mm	3/8	45 blue				RS156PA
22 mm	3/8	45 blue				RS194PA			
22 mm	3/8	90 blue					RS161PA		
24 mm	3/8	45 blue					RS204PA	RS121PA	
		75 blue					RS125PA	RS124PA	RS123PA

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle

length	needles curvature	thread length (cm) color	10/0	9/0	8/0	7/0	6/0	5/0	4/0
5 mm	1/2	13 black	RS231PA						
10 mm	1/2	75 blue				RS85PA			
12 mm	1/2	75 blue					RS83PA	RS84PA	
13 mm	1/2	75 blue						RS218PA	

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle

length	needles curvature		thread length (cm) color	6/0 5/0 4/0 3/0 2/0 0 1						
16 mm	1/2		75 blue		RS230PA	RS82PA	RS227PA			
16 mm	1/2		75 blue	RS171PA		RS132PA				
			90 blue		RS14PA	RS15PA	RS205PA			
17 mm	1/2		75 blue				RS193PA			
20 mm	1/2		75 blue		RS16PA	RS17PA	RS18PA	RS19PA	RS20PA	
20 mm	1/2		90 blue			RS105PA				
22 mm	1/2		75 blue			RS236PA	RS237PA			
			90 blue			RS167PA	RS180PA	RS150PA		RS149PA
26 mm	1/2		75 blue				RS24PA	RS25PA	RS97PA	RS127PA
			90 blue						RS26PA	
26 mm	1/2		90 blue					RS113PA		
			120 blue					RS104PA		
30 mm	1/2		75 blue		RS107PA	RS155PA		RS32PA	RS33PA	RS34PA
			90 blue					RS35PA	RS142PA	RS143PA
			100 blue					RS36PA	RS37PA	RS38PA
30 mm heavy	1/2		75 blue					RS102PA		
30 mm	1/2		150 blue					RS223PA	RS153PA	
35 mm	1/2		150 blue							RS137PA
36 mm heavy	1/2		45 blue						RS170PA	RS169PA

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament




round needle

length	needles curvature	thread length (cm) color	thread							
			5/0	4/0	3/0	2/0	0	1	2	
37 mm	1/2	75 blue				RS53PA	RS54PA	RS55PA	RS128PA	
		90 blue							RS56PA	
		100 blue						RS57PA	RS76PA	
40 mm	1/2	75 blue						RS60PA	RS61PA	RS62PA
		90 blue								RS110PA
		100 blue						RS63PA	RS64PA	RS65PA
40 mm heavy	1/2	75 blue								RS221PA
		100 blue						RS202PA		
40 mm	1/2	150 blue						RS139PA		
		200 blue						RS238PA		
40 mm heavy	1/2	150 blue						RS203PA		
48 mm	1/2	90 blue								RS108PA
		100 blue						RS189PA	RS66PA	RS67PA
		200 blue							RS165PA	
		240 blue							RS186PA	
48 mm heavy	1/2	100 blue						RS247PA		
48 mm	1/2	150 blue						RS190PA	RS181PA	
		240 blue							RS188PA	
48 mm heavy	1/2	150 blue						RS88PA		RS103PA
		200 blue							RS220PA	

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament




round needle

length	needles curvature	thread length (cm) color								
			5/0	4/0	3/0	2/0	0	1	2	
50 mm	1/2	100 blue							RS68PA	RS234PA
		150 blue						RS69PA	RS70PA	
50 mm heavy	1/2	100 blue							RS196PA	
50 mm	1/2	150 blue							RS233PA	
		240 blue loop							RS114PA	
60 mm	1/2	75 blue 2 vinyl tubes 60 mm						RS242PA		
		90 blue							RS241PA	
65 mm	1/2	150 blue							RS71PA	RS72PA
65 mm	1/2	120 blue						RS222PA		
		150 blue								RS192PA
		240 blue loop						RS187PA		
76 mm	1/2	75 blue						RS100PA	RS73PA	
80 mm	1/2	100 blue							RS166PA	

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



round needle

length	needles curvature	thread length (cm) color								
			5/0	4/0	3/0	2/0	0	1	2	
60 mm	straight	75 blue				RS148PA	RS91PA			
		90 blue					RS129PA	RS154PA		
60 mm	straight	90 blue					RS249PA			
65 mm	straight	90 blue					RS229PA			
70 mm	straight	75 blue					RS199PA			
70 mm	straight	75 blue					RS98PA			

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



cutting needle

length	needles curvature		thread length (cm) color	6/0 5/0 4/0 3/0 2/0 0 1							
8 mm cosmetic	3/8		45 blue	DK130PA							
12 mm	3/8		45 blue	DK71PA	DK95PA	DK67PA	DK85PA	DK86PA			
			75 blue	DK170PA	DK127PA	DK68PA					
12 mm cosmetic	3/8		75 blue		DK173PA						
16 mm	3/8		45 blue	DK01PA	DK02PA	DK03PA	DK87PA	DK88PA			
			45 black	DK174PA							
			75 blue	DK04PA	DK05PA	DK06PA	DK128PA				
16 mm cosmetic	3/8		45 blue	DK07PA	DK08PA	DK168PA					
			75 blue	DK09PA	DK10PA						
19 mm	3/8		45 blue	DK123PA	DK11PA	DK12PA	DK13PA	DK89PA			
			75 blue		DK14PA	DK15PA	DK16PA	DK129PA			
			90 blue		DK17PA	DK18PA	DK19PA	DK20PA			
19 mm cosmetic	3/8		45 blue		DK132PA	DK107PA	DK92PA				
24 mm	3/8		45 blue		DK100PA	DK21PA	DK22PA	DK23PA			
			75 blue			DK24PA	DK25PA	DK26PA	DK79PA		
			90 blue			DK27PA	DK28PA	DK29PA	DK172PA		
24 mm cosmetic	3/8		75 blue			DK154PA	DK146PA				
26 mm	3/8		45 blue			DK30PA	DK31PA	DK32PA			
			75 blue			DK33PA	DK34PA	DK35PA	DK69PA	DK165PA	
			90 blue			DK36PA	DK37PA	DK38PA	DK115PA		
26 mm cosmetic	3/8		45 blue				DK156PA				
			75 blue				DK114PA	DK112PA			

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.












NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



cutting needle

length	needles curvature	thread length (cm) color							
			5/0	4/0	3/0	2/0	0	1	2
30 mm	3/8 	45 blue			DK135PA	DK39PA	DK40PA		
		75 blue			DK41PA	DK42PA	DK43PA	DK75PA	
		90 blue			DK72PA	DK73PA	DK74PA	DK145PA	
35 mm	3/8 	75 blue			DK133PA	DK103PA			
		90 blue			DK163PA	DK84PA	DK105PA		
36 mm cosmetic	3/8 	75 blue			DK113PA				
39 mm	3/8 	45 blue						DK158PA	
		75 blue		DK93PA	DK94PA	DK78PA	DK46PA	DK47PA	
		90 blue				DK122PA	DK137PA	DK142PA	DK157PA
40 mm	3/8 	45 blue				DK48PA	DK49PA		
		75 blue				DK50PA	DK51PA	DK101PA	
		90 blue				DK150PA			
40 mm heavy	3/8 	45 blue				DK121PA			
		75 blue				DK125PA			
48 mm	3/8 	100 blue				DK140PA			
60 mm	3/8 	100 blue				DK99PA	DK76PA		DK77PA
90 mm	3/8 	45 blue						DK124PA	
		75 blue					DK153PA		
		90 blue						DK59PA	DK60PA
		100 blue				DK61PA	DK62PA	DK63PA	DK64PA

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.










NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



cutting needle

length	needles curvature	thread length (cm) color	thread						
			5/0	4/0	3/0	2/0	0	1	2
90 mm	3/8 	2x50 blue							DK155PA
90 mm	3/8 	75 blue							DK166PA
20 mm	1/2 	45 blue			DK181PA	DK182PA			
39 mm	1/2 	75 blue				DK136PA	DK97PA		
		90 blue					DK44PA	DK45PA	
40 mm	1/2 	45 blue					DK111PA		
		75 blue					DK110PA		
		90 blue						DK187PA	
		100 blue					DK54PA	DK180PA	
		150 blue						DK188PA	
48 mm	1/2 	75 blue						DK57PA	
		100 blue						DK58PA	
90 mm	1/2 	100 blue							DK65PA

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



cutting needle

length	needles curvature	thread length (cm) color								
			5/0	4/0	3/0	2/0	0	1	2	
51 mm	straight	75 blue				DG030PA	DG031PA	DG032PA		
		90 blue					DG040PA			
55 mm	straight	75 blue					DG07PA			
		100 blue					DG08PA			
60 mm	straight	75 blue				DG020PA	DG018PA	DG029PA		
		90 blue					DG021PA	DG041PA		
60 mm	straight	90 blue						DG26PA		
75 mm	straight	75 blue					DG019PA	DG043PA		

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature		thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
8 mm	3/8		45 blue	DKO01PA			DKO226PA	DKO244PA		
8 mm cosmetic	3/8		45 blue				DKO227PA			
			45 black				DKO379PA	DKO381PA		
11 mm	3/8		45 blue		DKO04PA	DKO05PA				
			45 black				DKO312PA			
			75 black					DKO206PA		
11 mm cosmetic	3/8		45 blue		DKO194PA	DKO235PA	DKO323PA			
12 mm	3/8		45 blue			DKO08PA	DKO09PA	DKO10PA	DKO327PA	
			75 blue			DKO11PA	DKO12PA	DKO13PA		
12 mm cosmetic	3/8		45 blue			DKO14PA	DKO15PA	DKO16PA		
			45 black			DKO318PA	DKO317PA			
			75 blue			DKO17PA	DKO18PA	DKO19PA		
13 mm	3/8		45 blue			DKO20PA	DKO21PA	DKO224PA		
			75 blue				DKO278PA			
13 mm	3/8		45 blue					DKO346PA		
13 mm cosmetic	3/8		45 blue			DKO213PA	DKO22PA	DKO251PA		
16 mm	3/8		45 blue			DKO23PA	DKO24PA	DKO25PA	DKO26PA	
			45 black				DKO361PA			
			75 blue			DKO27PA	DKO28PA	DKO29PA	DKO30PA	DKO220PA
16 mm cosmetic	3/8		45 blue			DKO31PA	DKO32PA	DKO33PA	DKO333PA	
			75 blue			DKO34PA	DKO35PA	DKO36PA		

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature	thread length (cm) color	6/0 5/0 4/0 3/0 2/0 0 1						
			6/0	5/0	4/0	3/0	2/0	0	1
19 mm	3/8	45 blue	DKO38PA	DKO39PA	DKO40PA	DKO41PA	DKO42PA	DKO358PA	
		45 black	DKO376PA	DKO360PA					
		75 blue		DKO43PA	DKO44PA	DKO45PA	DKO46PA	DKO331PA	
		75 black		DKO377PA					
		90 blue			DKO47PA	DKO48PA	DKO49PA		
19 mm cosmetic	3/8	45 blue		DKO50PA	DKO51PA	DKO52PA	DKO277PA		
		45 black		DKO316PA					
		75 blue	DKO53PA	DKO54PA	DKO55PA	DKO56PA			
22 mm	3/8	45 blue		DKO279PA	DKO321PA				
		75 blue			DKO276PA	DKO275PA	DKO263PA		
		90 blue				DKO192PA			
24 mm	3/8	45 blue	DKO234PA	DKO59PA	DKO60PA	DKO61PA	DKO62PA	DKO63PA	DKO301PA
		75 blue		DKO64PA	DKO65PA	DKO66PA	DKO67PA	DKO68PA	
		90 blue			DKO69PA	DKO70PA	DKO71PA		
24 mm cosmetic	3/8	45 blue		DKO236PA		DKO72PA	DKO73PA		
		75 blue				DKO74PA	DKO75PA		
26 mm	3/8	45 blue		DKO243PA	DKO189PA	DKO77PA	DKO78PA	DKO79PA	
		75 blue			DKO80PA	DKO81PA	DKO82PA	DKO83PA	DKO84PA
		90 blue			DKO285PA		DKO320PA	DKO292PA	
		100 blue							DKO85PA
26 mm cosmetic	3/8	45 blue			DKO382PA	DKO201PA	DKO294PA		
		75 blue			DKO282PA	DKO199PA	DKO200PA		

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature	thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
			30 mm	3/8	45 blue DKO373PA	DKO86PA	DKO87PA	DKO88PA	DKO190PA
		75 blue	DKO90PA	DKO91PA	DKO92PA	DKO93PA	DKO94PA	DKO215PA	
		90 blue	DKO256PA	DKO95PA	DKO96PA	DKO97PA	DKO98PA		
35 mm	3/8	45 blue		DKO105PA					
		75 blue		DKO106PA	DKO107PA	DKO205PA	DKO265PA		
		90 blue		DKO108PA	DKO109PA	DKO203PA			
37 mm	3/8	45 blue		DKO110PA					
		75 blue		DKO289PA	DKO219PA	DKO111PA	DKO112PA		
		90 blue		DKO370PA			DKO113PA		
39 mm	3/8	45 blue		DKO126PA	DKO208PA				
		75 blue		DKO127PA	DKO128PA	DKO129PA	DKO130PA	DKO369PA	
		90 blue		DKO131PA	DKO132PA	DKO133PA	DKO134PA	DKO135PA	
		100 blue				DKO286PA			
40 mm	3/8	45 blue				DKO204PA	DKO139PA		
		75 blue		DKO290PA	DKO140PA	DKO141PA	DKO142PA	DKO143PA	
		100 blue			DKO144PA	DKO145PA	DKO271PA	DKO191PA	
40 mm heavy	3/8	75 blue		DKO384PA	DKO341PA			DKO303PA	
40 mm cosmetic	3/8	75 blue				DKO221PA			
45 mm	3/8	75 blue				DKO157PA		DKO268PA	
		90 blue				DKO158PA	DKO306PA		
		100 blue				DKO160PA	DKO267PA		
45 mm cosmetic	3/8	75 blue					DKO385PA		

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature	thread length (cm) color								
			5/0	4/0	3/0	2/0	0	1	2	
48 mm	3/8	90 blue					DKO310PA		DKO272PA	DKO240PA
		100 blue					DKO319PA		DKO283PA	
60 mm	3/8	100 blue					DKO187PA	DKO171PA	DKO172PA	DKO173PA
77 mm	3/8	45 blue								DKO174PA
		75 blue							DKO264PA	DKO175PA
		90 blue							DKO375PA	
77 mm	3/8	75 blue							DKO238PA	DKO176PA
90 mm	3/8	2x50 blue							DKO326PA	
		75 blue								DKO269PA
		90 blue						DKO178PA		DKO179PA
		100 blue					DKO180PA	DKO181PA	DKO182PA	DKO183PA

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature		thread length (cm) color	thread								
				8/0	7/0	6/0	5/0	4/0	3/0	2/0		
8 mm	1/2		45 blue	DKO02PA				DKO228PA				
8 mm cosmetic	1/2		45 blue	DKO295PA								
16 mm	1/2		45 blue				DKO305PA	DKO366PA	DKO362PA	DKO367PA		
19 mm	1/2		45 blue						DKO262PA			
			90 blue						DKO186PA	DKO239PA		
22 mm	1/2		75 blue							DKO354PA	DKO355PA	
			90 blue							DKO58PA		
24 mm	1/2		75 blue							DKO237PA	DKO241PA	
			90 blue							DKO76PA		
25 mm	1/2		75 blue						DKO325PA			

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature	thread length (cm) color	thread								
			5/0	4/0	3/0	2/0	0	1	2		
26 mm	1/2	75 blue				DKO252PA	DKO254PA			DKO247PA	
30 mm	1/2	75 blue				DKO260PA	DKO99PA	DKO100PA	DKO101PA	DKO210PA	
		90 blue						DKO343PA			
37 mm	1/2	75 blue					DKO115PA	DKO116PA	DKO117PA	DKO249PA	
		90 blue					DKO118PA	DKO222PA	DKO119PA		
		100 blue							DKO120PA		
40 mm	1/2	75 blue						DKO149PA	DKO150PA	DKO151PA	
		90 blue						DKO233PA	DKO152PA	DKO153PA	
		100 blue						DKO154PA	DKO155PA	DKO156PA	
40 mm heavy	1/2	90 blue						DKO300PA			
48 mm	1/2	75 blue								DKO162PA	
		90 blue								DKO163PA	
		100 blue					DKO364PA	DKO322PA	DKO164PA	DKO242PA	
50 mm	1/2	75 blue								DKO165PA	
		90 blue								DKO372PA	DKO332PA
		100 blue								DKO166PA	
60 mm	1/2	100 blue					DKO253PA	DKO344PA			

notes



DKO

NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature	thread length (cm) color	thread							
			5/0	4/0	3/0	2/0	0	1	2	
45 mm	straight	75 blue			DG045PA					
51 mm	straight	90 blue		DG02PA	DG03PA	DG04PA	DG05PA	DG06PA		
55 mm	straight	100 blue				DG09PA				
60 mm	straight	75 blue		DG033PA	DG010PA	DG011PA				
		90 blue		DG044PA	DG012PA	DG013PA	DG014PA	DG015PA		
		100 blue				DG016PA	DG017PA			
60 mm	straight	75 blue		DG24PA	DG23PA	DG22PA				
		90 blue					DG042PA			
75 mm	straight	75 blue				DG039PA	DG028PA			

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



**cosmetic reverse cutting needle
with a precise blade (micro point)**

length	needles curvature		thread length (cm) color	6/0 5/0 4/0 3/0 2/0 0 1						
				6/0	5/0	4/0	3/0	2/0	0	1
13 mm	3/8		45 blue	ES23PA	ES13PA	ES29PA				
16 mm	3/8		45 blue	ES18PA	ES05PA	ES06PA	ES24PA			
			75 blue		ES19PA	ES25PA				
19 mm	3/8		45 blue		ES07PA	ES01PA	ES22PA			
			45 black		ES04PA					
			75 blue			ES08PA	ES09PA			
24 mm	3/8		75 blue			ES21PA	ES10PA	ES11PA		
26 mm	3/8		45 blue				ES03PA			
			75 blue				ES12PA	ES02PA		
30 mm	3/8		75 blue				ES32PA	ES14PA	ES27PA	
40 mm	3/8		75 blue					ES26PA		

notes



NYLON (polyamide)


surgical sutures, non-absorbable, monofilament



round needle with a cutting edge (sharpened)

length	needles curvature	thread length (cm)	color	5/0	4/0	3/0	2/0	0	1	2
24 mm	3/8	75	blue		DS16PA	DS17PA	DS18PA			
26 mm	3/8	75	blue				DS09PA			
16 mm	1/2	90	blue	DS10PA						
26 mm	1/2	75	blue				DS13PA			
26 mm	1/2	90	blue			DS11PA				
30 mm	1/2	75	blue						DS21PA	
		90	blue				DS05PA		DS14PA	
40 mm	1/2	75	blue						DS07PA	DS22PA
		100	blue					DS06PA		
45 mm	1/2	100	blue					DS08PA		
		200	blue						DS04PA	
45 mm	1/2	200	blue						DS20PA	
45 mm heavy	1/2	200	blue							DS23PA

notes



NYLON (polyamide)

surgical sutures, non-absorbable, monofilament



reverse cutting P needle blade

length	needles curvature	thread length (cm) color								
			6/0	5/0	4/0	3/0	2/0	0	1	
13 mm	3/8	75 blue			DM12PA					
16 mm	3/8	45 blue			DM10PA	DM01PA				
		45 black			DM09PA					
19 mm	3/8	75 blue			DM11PA	DM02PA	DM03PA			
		45 blue				DM04PA				
24 mm	3/8	75 blue				DM05PA	DM06PA			
		45 blue					DM08PA			
		75 blue				DM07PA				

notes

NYLON (polyamide)

surgical sutures, non-absorbable, monofilament

LIGATURES

threads length	color thread	7/0	6/0	5/0	4/0	3/0	2/0	0	1	2
2 x 75 cm	blue						P05PA		P06PA	
2 x 100 cm	blue					P011PA	P038PA	P012PA	P013PA	P014PA
3 x 45 cm	blue					P59PA	P49PA	P60PA		
6 x 45 cm	blue					P031PA	P032PA	P25PA	P028PA	P026PA
75 cm	blue								P03PA	P04PA
100 cm	blue						P07PA	P08PA	P09PA	P010PA
150 cm	blue					P016PA	P020PA	P017PA	P018PA	P019PA
150 cm in reel	blue							Pr017PA	Pr018PA	Pr019PA
180 cm	blue						P46PA	P45PA	P43PA	P44PA

notes