



# VISCOJECT™ **easy**

# VISCOJECT™ **eco**

INJECTION SYSTEMS

More for less



Medicel AG

Dornierstrasse 11

CH-9423 Altenrhein

SWITZERLAND

Tel. +41 71 727 10 50

Fax +41 71 727 10 55

[info@medicel.com](mailto:info@medicel.com)

[www.medicel.com](http://www.medicel.com)

**medicel**  
SWISS TECHNOLOGY FOR SURGERY



# VISCOJECT™ **easy**

VISCOJECT™ **easy** is the latest member of the VISCOJECT™ family. It simplifies the loading of the lens with its user-friendly guiding rails within the loading chamber.

The loading of the lens is performed with the aid of a loading unit, which defines the angular position of the cartridge wing and considerably simplifies the loading process.

The design of the cartridge tip is equivalent to the advanced VISCOJECT™-BIO cartridge.

## Single-use VISCOJECT™ **easy** injection system



- For hydrophilic and hydrophobic 1-piece lenses
- Simple handling with the aid of a loading unit
- Easy loading due to guiding rails within cartridge
- With classic gliding agent
- Same cartridge design as VISCOJECT™-BIO
- Incision sizes range from 2.9mm to 2.0mm

## Reusable VISCOJECT™ **easy** injection system



- High grade titanium
- Available as syringe-type or screw-type injector
- To be used with fully single-use cartridge sets



## VISCOJECT™ **eco**

Are you located in a price sensitive market?

The attractive price-performance ratio of the VISCOJECT™ **eco** is the perfect solution. Most importantly, equipped with the traditional VISCOJECT™ cartridge, the VISCOJECT™ **eco** does not compromise on quality and safety.

### Single-use VISCOJECT™ **eco** injection system



- For hydrophilic and hydrophobic 1-piece lenses
- Complete control of every step of the lens injection
- With soft tip technology
- With classic gliding agent
- Incision sizes range from 2.7mm to 2.4mm



- Cost-efficient packaging without compromising product quality and safety

single-use VISCOJECT™ **easy** Injector Set

| ITEM CODE | PRODUCT NAME                  | INCISION SIZE (mm) |                |              | IOL TYPE       |
|-----------|-------------------------------|--------------------|----------------|--------------|----------------|
|           |                               | wound-assisted     | into-the-wound | into-the-bag |                |
| VY2700    | VISCOJECT™ <b>easy</b> 2.7-1P | 2.2 - 2.3          | 2.5 - 2.6      | 2.9          | 1-piece lenses |
| VY2400    | VISCOJECT™ <b>easy</b> 2.4-1P | 2.1                | 2.4            | 2.7          | 1-piece lenses |
| VY2200    | VISCOJECT™ <b>easy</b> 2.2-1P | 1.9 - 2.0          | 2.2            | 2.5          | 1-piece lenses |

reusable VISCOJECT™ **easy** Injection System

| ITEM CODE | PRODUCT NAME                        | TYPE         | TO BE USED WITH VISCOJECT™ <b>easy</b> CARTRIDGE SETS |
|-----------|-------------------------------------|--------------|---|
| LI604305  | VISCOJECT™ <b>easy</b> „ATRAUMATIC“ | syringe type | VC2700 / VC2200                                       |
| LI604315  | VISCOJECT™ <b>easy</b> „ATRAUMATIC“ | screw type   | VC2700 / VC2200                                       |

single-use VISCOJECT™ **easy** Cartridge Set

| ITEM CODE | PRODUCT NAME                  | INCISION SIZE (mm) |                |              | IOL TYPE       |
|-----------|-------------------------------|--------------------|----------------|--------------|----------------|
|           |                               | wound-assisted     | into-the-wound | into-the-bag |                |
| VC2700    | VISCOJECT™ <b>easy</b> 2.7-1P | 2.2 - 2.3          | 2.5 - 2.6      | 2.9          | 1-piece lenses |
| VC2400    | VISCOJECT™ <b>easy</b> 2.4-1P | 2.1                | 2.4            | 2.7          | 1-piece lenses |
| VC2200    | VISCOJECT™ <b>easy</b> 2.2-1P | 1.9 - 2.0          | 2.2            | 2.5          | 1-piece lenses |

single-use VISCOJECT™ **eco** Injector Set

| ITEM CODE | PRODUCT NAME              | INCISION SIZE (mm) |                |              | IOL TYPE       |
|-----------|---------------------------|--------------------|----------------|--------------|----------------|
|           |                           | wound-assisted     | into-the-wound | into-the-bag |                |
| VE2700    | VISCOJECT™ <b>eco</b> 2.7 | 2.3                | 2.6 - 2.7      | 2.9 - 3.0    | 1-piece lenses |
| VE2200    | VISCOJECT™ <b>eco</b> 2.2 | 2.1 - 2.2          | 2.4            | 2.7 - 2.8    | 1-piece lenses |