


Orizon

Preloaded


Hydrophobic Acrylic
Foldable Intraocular Lens




Orizon

Preloaded

Hydrophobic Acrylic
Foldable Intraocular Lens



SWISS QUALITY
meditel
SWISS TECHNOLOGY FOR SURGERY



Premium Aspheric Optic Design that provides sharper vision


360° Square Edge surface that reduces PCO rates

Superb quality of raw materials with extremely low inflammatory response and optimal biocompatibility

Product Features


- Cast molding technology
- Superb quality of raw materials using pre-molded blanks
- Biocompatible hydrophobic acrylic material with UV absorption
- With Biconvex surfaces which minimize effects of tilt, decentration and aberration
- Preloaded with premium quality hydrophobic aspheric acrylic C Loop 6X13mm foldable intraocular lens.
- Notches for smooth manipulation
- Ergonomic grip for enhanced stability
- Absolute clarity throughout the delivery

SDHBP
Aspheric,
360° Square Edge



Optic Design: Biconvex
Optic Size: 6.00 mm
Overall Length: 13.00 mm *
Haptic Angle: 0°
Refractive Index : 1.56
A.Constant : 118.5

SDHBPY 
Aspheric,
360° Square Edge
Yellow, Blue Light Filtering



Optic Design: Biconvex
Optic Size: 6.00 mm
Overall Length: 13.00 mm *
Haptic Angle: 0°
Refractive Index : 1.49
A.Constant : 118.5

SIDAPHARM

21, Stageiriti & 24, Em. Fili Str.
Thessaloniki, GR- 54352, Greece, European Union
T +30 2310:906660 info@sidapharm.gr
F: +30 2310 989846 www.sidapharm.gr 

* additional sizes also available