

Ophthalmic Catalogue

Table of contents

Who we are 1

Intraocular Lenses 2 - 5

Microsurgical Ophthalmic Knives 6 - 7

Ophthalmic Cannulas 8

Ophthalmic Viscoelastics (OVDs) 9 - 10

Ophthalmic Solutions 10 - 11

Miscellaneous Disposables 12

Cataract Kits 13

Surgical Apparel 13

Who we are

Sidapharm was founded in 2000 and is a manufacturer of a vast range of CE-certified ophthalmic products such as intraocular lenses, knives and cannulas, viscoelastics, BSS, cataract kits and ophthalmic instruments. Based in Thessaloniki, Greece, a city which connects Europe with Asia, Middle East and Africa, the company managed to take advantage of its location, expand its network and build strong collaborations through these continents.

Nowadays, Sidapharm has presence in more than 100 countries, mostly through very dynamic and reliable distributors. Our focus is on investing in long-lasting partnerships and assisting our distributors in every possible way to become more competitive and capture a bigger market share.




Intraocular Lenses

Preloaded

Hydrophobic, Aspheric, One-Piece,
Acrylic, Foldable, Intraocular Lens

Orizon 
PRELOADED





PREMIUM PRECISION


- Cast molding technology
- Superb quality of raw materials using pre-molded blanks
- Biocompatible hydrophobic acrylic material with UV absorption
- With Biconvex surfaces which minimize effects of tilt, decentration and aberration
- 360°square-edge profile for reduced PCO rates
- Notches for smooth manipulation

Hydrophobic

One-Piece, Acrylic,
Foldable, Intraocular Lens



*raw material by
Contamac
EMPOWERING YOUR VISION

<p>SDHB6130 Aspheric, * 360° Square Edge</p> <p>*Also available in Spherical Design SDHB6130S</p>		<p>Optic Design: Biconvex Optic Size: 6.00 mm Overall Length: 13.00 mm</p>
<p>SDHB Aspheric, 360° Square Edge</p>		<p>Optic Design: Biconvex Optic Size: 6.00 mm Overall Length: 12.50 mm</p>
<p>SDHBY  Aspheric, 360° Square Edge Yellow, Blue Light Filtering</p>		<p>Optic Design: Biconvex Optic Size: 6.00 mm Overall Length: 12.50 mm</p>

Intraocular Lenses

Hydrophilic

One-Piece, Acrylic,
Foldable, Intraocular Lens




- Exceptional biocompatibility and UV absorption
- High Resolution, Diopter Accuracy
- Easy lens manipulation
- Controlled surgical handling

raw material by
Contamac
EMPOWERING YOUR VISION

SDA


Aspheric,
360° Square Edge,



Optic Design: Biconvex
Optic Size: 6.00 mm
Overall Length: 12.50 mm

SDAC


Aspheric,
360° Square Edge



Optic Design: Biconvex
Optic Size: 6.00 mm
Overall Length: 12.50 mm

SDACY


Aspheric,
360° Square Edge
Yellow, Blue Light Filtering



Optic Design: Biconvex
Optic Size: 6.00 mm
Overall Length: 12.50 mm

SDA4

Aspheric,
360° Square Edge



Optic Design: Biconvex
Optic Size: 6.00 mm
Overall Length: 11.00 mm

PMMA



Clinical quality PMMA with UV absorber


YAG laser compatible

Proven Biocompatibility

ANTERIOR CHAMBER

SDANT

4 points of fixation




Optic Design: Equiconvex
Optic Size: 6.00 mm
Overall Length: 12.50 mm
Haptic Design: S Loop

POSTERIOR CHAMBER

SDALP

2 Dialing holes



Optic Design: Equiconvex
Optic Size: 6.00 mm
Overall Length: 12.50 mm
Haptic Design: C Loop

* PMMA SDHS model produced on demand

Ophthalmic Microsurgical Knives


- Electropolishing technology
- Minimized micro-roughness
- Cleanability of surfaces
- Low contamination levels

Effects of electropolishing on the surface finish

Surface Before Electropolishing

Surface After Electropolishing

Knife Comparison



Precision

Diamond like sharpness
Exact Dimensions


Performance

Incisions with accuracy

Quality

Matte anti-glare

Sharpness Comparison



Tip penetration chart, exemplifies exceptional penetration during procedure, using minimal force

Ophthalmic Microsurgical Knives



Premium quality knives, deliver excellent levels of sharpness, precision and enhanced control while maintaining low cost.



Stab Knives



Ref. No.	Description
62000	15°
62021	22.5°
62001	30°
62002	45°

Crescent Knives



Ref. No.	Description
62003	Angled
62003-1	Straight

Slit Knives



Ref. No.	Description
62016	2.20mm Angled
62017	2.40mm Angled
62052	2.60mm Angled
62018-1	2.75mm Angled
62018	2.80mm Angled
62019	3.00mm Angled
62020	3.20mm Angled

Slit Knives - Double Bevel



Ref. No.	Description
62016-DB	2.20mm Angled
62017-DB	2.40mm Angled
620181-DB	2.75mm Angled
62018-DB	2.80mm Angled
62019-DB	3.00mm Angled
62020-DB	3.20mm Angled

MVR Knives




Ref. No.	Description
62004	19G Straight
62009	19G Angled
62005	20G Straight
62010	20G Angled
62007	23G Straight
62012	23G Angled

Ophthalmic Cannulas

- Steel and Plastic of the highest quality
- Guaranteed hygiene
- Constant Development & Improvement

CYSTOTOMES


Formed, 16mm tip



Ref. No.	Description
77006	25G, 16mm
77007	27G, 16mm

ANTERIOR CHAMBER

Viscoelastic, angled 22mm, formed



77104*	25G, 22mm
77082	27G, 22mm

Rycroft / Air Cannula , angled tip



77044*	23G, 4mm
77001	25G, 8mm
77003 - 1*	27G, 4mm
77056*	27G, 8mm
77003	27G, 10mm
77046*	30G, 4mm
77004	30G, 6mm

HYDRODISSECTORS

Angled, flat tip



77062*	23G, 8mm tip
77010	25G, 8mm tip
77012	27G, 8mm tip
77013*	30G, 8mm tip

Curved, flat tip



77014	27G, Horizontal tip
77063*	27G, Vertical tip

ANESTHESIA

Retrobulbar, straight Atkinson tip



Ref. No.	Description
77027	25G, 38mm tip

Subtenon




77032	19G, 22mm tip
-------	---------------

Peribulbar, straight Atkinson tip




77030*	25G, 22mm tip
77038*	25G, 32mm tip

A/C MAINTAINERS



77015*	19G, Screw tip
77016*	20G, Screw tip

SIMCOE



77022*	20G
77025	23G

LACRIMAL



77021	26G	28mm
-------	-----	------

* produced on demand

Ophthalmic Viscoelastics (OVDs)

Sodium Hyaluronate

Ophthalmic Solutions

Product Features :

- Easy to inject and remove from the eye
- Latex-free
- Provides excellent endothelial protection
- Ideal for cartridge lubrication during IOL implantation
- No refrigeration required



Reference Numbers per concentration and volume


	1.0 %	1.4 %	1.6 %	1.8 %	2.0 %	3.0 %
1.0 mL	10000	10005	10006	10002	10009	10007
1.5 mL*	10010	10014	10016	10018	-	10030

Available in the following concentrations


Viscoelastics provide the protection and control needed during ophthalmic procedures of :

- Cataract extraction
- IOL implantation
- Glaucoma filtering surgery
- Corneal transplantation surgery


SYisc 1.0 %




SYisc 1.4 %



SYisc 1.6 %



SYisc 1.8 %



SYisc 2.0 %



SYisc 3.0 %



*available upon request

Ophthalmic Viscoelastics (OVDs)


SHPMC 2%

Offering an ophthalmic viscosurgical device that has low viscosity properties and is ideal for coating and protecting the eye during anterior segment ophthalmic procedures.

HPMC 2% exhibits moderate viscoelasticity characteristics, with retentive but also excellent tissue coating abilities. Offers good space maintenance and excellent tissue protection throughout the ophthalmic procedure.

KEY Product features :

- Low molecular weight and low surface tension
- High concentration of HPMC
- Non-antigenic (allergic)
- Non-pyrogenic
- Completely transparent
- Sterile pre-filled syringe supplied with a cannula
- Available in 2 ml luer-lock glass syringe
- Preservative-free
- Ambient Storage



Ref. No. : 10003

Ophthalmic Solutions

BSS Balanced Salt Solution

BSS is a sterile physiological solution for use in irrigating ocular tissues of the eye. It hydrates and protects tissues and cells during any ophthalmic operation.

Suitable for :

- Refractive surgery
- Phacoemulsification
- Vitrectomy
- Corneal irrigation

High quality non collapsible plastic

Available in a two way volumetric 500ml plastic bottle with retractable, attached, hanging loop at the base.



Ref. No. : 76000

Trypan Blue

Ophthalmic Solution 0.06%

- Excellent visualization of the capsulorhexis in eyes with mature cataract or small pupils
- Clear outline of the capsulorhexis rim during the entire surgery
- Capsule staining is quick, easy and safe to perform
- Creates excellent staining contrast for anterior and posterior capsulorhexis
- Reduces surgical time

SBlue PFS - 1ml Prefilled syringe

- Ready to use syringe
- Ergonomic grip for better control of the injection
- Disposable ophthalmic cannula included



Ref. No. : 84001

SBlue - 1ml vial



Ref. No. : 84000

INDICATED FOR :

- Cataract surgery with poor red reflex
- Phacoemulsification
- Small incision cataract surgery
- Extra capsular cataract surgery
- Incomplete capsulorhexis
- Anterior/posterior capsulotomy
- Any high risk cataract case
- Trauma Cases
- Teaching trainee surgeons

- Easy application
- Reduces duration of surgery

Miscellaneous Disposables

Fluorescein Sodium

Ophthalmic Strips
1.0 mg


Uniform distribution ensuring optimal staining

Enhanced detection of :

- Corneal abrasions
- Foreign bodies
- Herpetic corneal ulcers

Individually packed, sterile ophthalmic strips

Ref. No. : 91000




Eye Spears

A soft, absorbing tool to drain excessive fluids

Provided in a sterile pack of 5 pieces

PVA material



Product Features :

Lint free

High density formulation

Rapid fluid absorption

Gamma radiation sterilized





Double pouch packing

Ref. No. :70001



Surgical Apparel

After years of constant research and cooperation with doctors our complete line of surgical apparel is here to cover your needs. Individual sterile packaging optimizes time and resources.

	<p>20014 Ophthalmic Drape 45g SMS With pouch and incision film 100cm x 120 cm</p>		<p>20028 Ophthalmic Drape With aperture 50cm x 60cm</p>
	<p>20039 Ophthalmic Drape With pouch and incision film 50cm x 60 cm</p>		<p>25001 Surgical Gown 50g SMS non-woven Fluid repellent Provides unrestricted movement</p>

Cataract Kits


Specially designed for cataract operations
Ideal as an one-stop supply solution
Excellent quality for great price

Basic Cataract Pack


Designed to be user-friendly, convenient and efficient in order to meet any special needs of the doctor and the surgical team.

Contents :

- 1 x Table Cover (140cm x 150cm)
- 1 x Ophthalmic Drape with pouch (100cm x 120cm)
- 1 x Surgical Drape (100cm x 160cm)
- 1x 1ml syringe 26g
- 2x 3ml syringe
- 2x 5ml syringe
- 1x 10ml syringe
- 10x Gauze swab 7.5 x 7.5
- 1x Eye Pad
- 2x Gowns with towels (Large)
- 1x 60 ml plastic bowl
- 5x Eye spears





21, Stageiriti & 24, Em. Fili Str.
Thessaloniki, GR- 54352, Greece, European Union
T: +30 2310 906660 info@sidapharm.gr
F: +30 2310 989846 www.sidapharm.gr 

All figures and images used are for illustrative purposes only. They are indicative of the quality and style of the specification and may not represent wholly the actual items. The specifications may be subject to change. Please ask our sales adviser for details. Whilst all reasonable care has been taken in providing true information, we accept no responsibility for the accuracy of any information contained herein.

Ver 4.0 6/2020