



POLYPROPYLENE

POLYPROPYLENE

Non-absorbable, synthetic surgical monofilament thread



- Polypropylene
- Uncoated
- Thread may be blue or black
- Size 10/0 – 2 USP
- Features
 - Uniform structure with smooth surface
 - Tensile strength
- Indications
 - Fascia and skin surgery
 - General, pediatric, plastic surgery
 - Suturing and integration of blood vessels, vascular prostheses
 - Cardiac surgery
 - Neurosurgery
 - Ophthalmology
 - Surgical mesh implantation
- Period of validity 5 years

yavo

YAVO Sp. z o. o. 97-400 Bełchatów, ul. Bawełniana 17, Polska
tel.: +48 44 632 94 11, fax:+48 44 632 49 17,
e-mail: info@yavo.com.pl, www.yavo.com.pl



POLYPROPYLENE



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



spatula (lanced needle)

length	needles		thread length (cm) color	10/0	9/0	8/0	7/0	6/0	5/0	4/0
	curvature	diameter								
6 mm	3/8		15 black	LZ28PP						
			23 black		LZ09PP					
			30 black	LZ01PP						
6 mm micropoint	3/8	220 µm	23 black		LZ25PP					
6 mm	3/8		30 black	LZ02PP	LZ03PP					
			30 blue	LZ29PP						
6 mm micropoint	3/8		30 black	LZ30PP						
6 mm micropoint	3/8		30 blue	LZ45PP						
6,2 mm	3/8		8 black		LZ16PP					
6,2 mm	3/8		30 black	LZ15PP			LZ08PP	LZ07PP		
6,5 mm micropoint	3/8		23 black	LZ27PP						
9 mm	3/8		30 blue				LZ17PP			
5,5 mm	1/2		30 black	LZ13PP						
16 mm	straight		23 black	LZ06PP						
			25 black	LZ10PP						
16 mm micropoint	straight		23 black	LZ23PP						
			23 blue	LZ26PP						
			25 black	LZ14PP						
16 mm	straight	150 µm	23 black	LZ24PP						

notes




POLYPROPYLENE

surgical sutures, non-absorbable, monofilament



round needle

length	needles		thread length (cm) color	10/0	9/0	8/0	7/0	6/0	5/0	4/0
	curvature	diameter								
4 mm	3/8		15 black	RS345PP						
5 mm	3/8		15 black	RS266PP	RS419PP					
5 mm	3/8		15 blue	RS555PP						
5 mm	3/8		30 black	RS361PP						
6 mm	3/8		60 blue		RS554PP					
6 mm	3/8		45 blue		RS556PP					
			30 black		RS362PP					
			45 blue			RS265PP				
			60 blue			RS344PP	RS451PP			
6 mm	3/8		60 blue			RS445PP				
6,5 mm	3/8		13 blue			RS349PP				
8 mm	3/8	140 µm	15 blue			RS219PP				
			60 blue			RS03PP	RS209PP	RS299PP		
8 mm	3/8		60 blue			RS05PP	RS06PP	RS284PP	RS437PP	
			75 blue		RS375PP	RS07PP	RS300PP			
9 mm	3/8		75 blue				RS210PP	RS431PP		
9 mm	3/8		45 blue			RS08PP	RS09PP	RS10PP		
			60 blue				RS11PP	RS12PP		
			75 blue			RS13PP	RS250PP	RS14PP		
9,3 mm	3/8	220 µm	45 blue					RS418PP		
			60 blue			RS15PP	RS206PP	RS249PP		
			75 blue				RS294PP			
10 mm	3/8	200 µm	60 blue			RS26PP	RS325PP	RS27PP		
			75 blue				RS226PP	RS220PP		
11 mm	3/8		13 blue				RS408PP			
			60 blue					RS303PP		
			75 blue				RS399PP			



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



round needle

length	needles curvature		thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
12 mm	3/8		60 blue		RS252PP	RS31PP	RS32PP			
			75 blue			RS33PP	RS34PP	RS438PP		
12 mm	3/8		75 blue				RS279PP	RS396PP		
13 mm	3/8		45 blue		RS39PP	RS230PP				
			60 blue		RS308PP	RS307PP				
			75 blue		RS248PP	RS40PP	RS41PP	RS42PP		
			90 blue			RS247PP				
13 mm	3/8		60 blue		RS44PP	RS45PP	RS46PP			
			75 blue		RS47PP	RS48PP	RS49PP	RS50PP		
			90 blue			RS304PP	RS260PP			
17 mm	3/8		45 blue			RS261PP				
			60 blue					RS259PP		
			75 blue				RS90PP	RS91PP		
			90 blue				RS239PP	RS277PP		
17 mm	3/8		75 blue					RS390PP	RS371PP	
			90 blue					RS229PP		
20 mm	3/8		60 blue					RS435PP		
			75 blue					RS235PP		
22 mm	3/8		75 blue						RS372PP	RS373PP
22 mm	3/8		75 blue						RS441PP	
9 mm	1/2		75 blue		RS377PP					
9,3 mm	1/2		60 blue	RS16PP	RS17PP	RS18PP				

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



round needle

length	needles curvature		thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
10 mm	1/2		60 blue				RS245PP			
			75 blue	RS434PP	RS297PP	RS340PP	RS274PP			
10 mm	1/2		60 blue		RS290PP	RS244PP	RS243PP			
			75 blue		RS29PP	RS30PP				
13 mm	1/2		60 blue		RS51PP	RS52PP				
			75 blue		RS53PP	RS54PP	RS55PP	RS391PP		
13 mm	1/2		60 blue		RS56PP	RS269PP	RS285PP	RS57PP		
			75 blue			RS58PP	RS59PP	RS60PP		
			90 blue			RS61PP	RS62PP	RS379PP		
16 mm	1/2		45 blue				RS69PP			
			75 blue			RS70PP	RS71PP	RS72PP	RS73PP	
16 mm thin	1/2		75 blue				RS326PP	RS327PP	RS295PP	
16 mm	1/2		45 blue			RS242PP	RS77PP			
			75 blue			RS78PP	RS79PP	RS80PP	RS81PP	
			90 blue		RS268PP	RS82PP	RS83PP	RS84PP	RS85PP	
17 mm	1/2		75 blue			RS313PP	RS86PP	RS87PP	RS88PP	
			90 blue					RS332PP	RS236PP	
17 mm	1/2		60 blue				RS394PP			
			75 blue				RS92PP	RS93PP	RS283PP	
			90 blue			RS94PP	RS95PP	RS96PP	RS97PP	RS334PP
20 mm	1/2		45 blue					RS98PP	RS99PP	RS100PP
			75 blue				RS101PP	RS102PP	RS103PP	RS104PP
20 mm	1/2		75 blue				RS207PP	RS107PP	RS108PP	
			90 blue				RS289PP	RS109PP	RS110PP	RS111PP

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



POLYPROPYLENE

surgical sutures, non-absorbable, monofilament




round needle

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
22 mm	1/2		45 blue			RS115PP				
			75 blue		RS116PP	RS117PP	RS118PP	RS119PP	RS256PP	
22 mm heavy	1/2		75 blue			RS273PP				
22 mm	1/2		75 blue	RS122PP	RS123PP	RS124PP				
			90 blue		RS125PP	RS126PP	RS127PP			
25 mm	1/2		75 blue		RS350PP	RS129PP	RS130PP	RS439PP		
			90 blue		RS131PP	RS132PP	RS133PP			
25 mm	1/2		75 blue		RS134PP	RS135PP	RS136PP	RS137PP		
			90 blue		RS138PP	RS139PP	RS140PP	RS141PP		
26 mm thin	1/2		75 blue			RS328PP	RS275PP	RS296PP		
26 mm	1/2		75 blue		RS142PP	RS143PP	RS144PP	RS145PP	RS301PP	
			90 blue		RS356PP	RS315PP	RS146PP	RS147PP		
			120 blue			RS212PP	RS150PP			
26 mm	1/2		75 blue	RS363PP	RS151PP	RS152PP	RS153PP	RS154PP		
			90 blue	RS393PP	RS155PP	RS156PP	RS157PP	RS158PP		
30 mm	1/2		75 blue			RS161PP	RS162PP	RS163PP	RS164PP	
			90 blue			RS165PP	RS166PP	RS167PP		
			100 blue					RS168PP	RS436PP	
30 mm heavy	1/2		75 blue				RS272PP			
			100 blue					RS169PP	RS170PP	

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



round needle

length	needles curvature	thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
30 mm	1/2	75 blue				RS171PP			
		90 blue			RS172PP	RS173PP			
35 mm	1/2	75 blue				RS341PP	RS174PP	RS175PP	
		100 blue				RS311PP	RS177PP	RS178PP	RS440PP
37 mm	1/2	75 blue				RS179PP	RS180PP	RS181PP	RS232PP
		90 blue						RS291PP	
		100 blue						RS182PP	
37 mm	1/2	90 blue			RS286PP	RS214PP		RS443PP	
		120 blue				RS321PP			
37 mm heavy	1/2	75 blue				RS346PP	RS329PP	RS218PP	RS217PP
40 mm	1/2	90 blue					RS263PP	RS264PP	
40 mm	1/2	75 blue				RS228PP	RS257PP	RS187PP	RS267PP
		90 blue				RS287PP	RS188PP	RS224PP	
		100 blue					RS189PP	RS190PP	
40 mm heavy	1/2	75 blue						RS330PP	
		100 blue					RS221PP	RS254PP	RS395PP
45 mm	1/2	75 blue						RS191PP	
		100 blue					RS319PP		

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



POLYPROPYLENE

surgical sutures, non-absorbable, monofilament




round needle

length	needles curvature	thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
48 mm	1/2	75 blue						RS192PP	RS416PP
		90 blue						RS320PP	
		100 blue					RS193PP		RS194PP
		150 blue							RS368PP
48 mm heavy	1/2	90 blue						RS432PP	
48 mm	1/2	150 blue					RS339PP	RS430PP	
65 mm	1/2	150 blue						RS197PP	RS198PP
65 mm heavy	1/2	100 blue						RS444PP	
65 mm	1/2	240 blue					RS227PP	RS309PP	
30 mm	hook	100 blue				RS281PP			
51 mm	straight	90 blue				RS316PP			
51 mm	straight	90 blue				RS199PP			

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.





POLYPROPYLENE

surgical sutures, non-absorbable, monofilament



round needle

length	needles curvature	thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
60 mm	straight 	75 blue			RS200PP	RS201PP	RS446PP		
65 mm	straight 	75 blue				RS202PP			
70 mm	straight 	75 blue				RS205PP			
		100 blue				RS258PP			

notes



POLYPROPYLENE

surgical sutures, non-absorbable, monofilament



black needle

length	needles curvature		thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
9 mm	3/8		60 blue	RSC01PP	RSC02PP	RSC03PP				
			75 blue			RSC04PP				
9,5 mm	3/8		60 blue		RSC05PP					
13 mm	3/8		75 blue				RSC18PP			
13 mm	3/8		60 blue		RSC06PP	RSC07PP	RSC08PP			
			75 blue			RSC09PP				
			90 blue				RSC10PP			
16 mm	1/2		90 blue				RSC11PP			
17 mm	1/2		90 blue					RSC12PP		
20 mm	1/2		90 blue					RSC13PP		
22 mm	1/2		90 blue					RSC14PP		
25 mm	1/2		90 blue						RSC15PP	RSC16PP
26 mm	1/2		90 blue					RSC17PP		

notes



POLYPROPYLENE

surgical sutures, non-absorbable, monofilament



cutting needle

length	needles curvature		thread length (cm) color	6/0	5/0	4/0	3/0	2/0	0	1
11 mm cosmetic	3/8		45 blue	DK59PP	DK63PP					
13 mm cosmetic	3/8		45 blue			DK62PP				
16 mm	3/8		45 blue	DK01PP	DK02PP	DK03PP				
			75 blue	DK04PP	DK05PP	DK06PP				
16 mm cosmetic	3/8		45 blue		DK37PP	DK42PP				
19 mm	3/8		45 blue	DK30PP	DK07PP	DK08PP	DK29PP			
			75 blue	DK09PP	DK10PP	DK11PP	DK33PP			
24 mm	3/8		45 blue			DK12PP	DK13PP	DK35PP		
			75 blue			DK14PP	DK15PP	DK31PP		
24 mm cosmetic	3/8		75 blue					DK60PP		
26 mm	3/8		45 blue				DK16PP			
			75 blue				DK17PP		DK28PP	
26 mm cosmetic	3/8		45 blue				DK40PP			
30 mm	3/8		75 blue				DK26PP	DK27PP	DK36PP	
35 mm cosmetic	3/8		75 blue				DK50PP			
39 mm	3/8		75 blue				DK25PP	DK18PP	DK41PP	

notes




POLYPROPYLENE

surgical sutures, non-absorbable, monofilament



cutting needle

length	needles curvature	thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
90 mm	3/8 	90 blue							DK48PP
		100 blue					DK52PP	DK22PP	DK23PP
90 mm	3/8  with vinyl tube 40 mm	50 blue						DK49PP	
		2x50 blue						DK24PP	
26 mm heavy	1/2 	100 blue						DK64PP	
55 mm	straight 	90 blue				DG03PP			
60 mm	straight 	75 blue			DG04PP	DG05PP	DG09PP		
		90 blue				DG06PP			

notes



POLYPROPYLENE

surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature		thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
8 mm	3/8		45 blue	DKO119PP						
			75 blue	DKO146PP						
8 mm cosmetic	3/8		45 blue			DKO116PP				
11 mm	3/8		45 blue			DKO01PP	DKO93PP			
11 mm cosmetic	3/8		45 blue		DKO81PP	DKO130PP				
12 mm	3/8		45 blue			DKO02PP	DKO104PP		DKO179PP	
			75 blue		DKO148PP	DKO147PP	DKO161PP	DKO98PP		
12 mm cosmetic	3/8		45 blue			DKO03PP	DKO96PP			
16 mm	3/8		45 blue			DKO91PP	DKO07PP	DKO08PP		
			75 blue		DKO135PP	DKO155PP	DKO09PP	DKO10PP	DKO11PP	
16 mm cosmetic	3/8		45 blue			DKO107PP	DKO12PP	DKO13PP	DKO79PP	
19 mm	3/8		45 blue				DKO15PP	DKO16PP	DKO17PP	
			75 blue				DKO18PP	DKO19PP	DKO20PP	
19 mm cosmetic	3/8		45 blue				DKO82PP	DKO21PP	DKO22PP	
			75 blue					DKO83PP		

notes



POLYPROPYLENE

surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
24 mm	3/8		45 blue	DKO177PP	DKO26PP	DKO27PP	DKO28PP	DKO29PP		
			75 blue		DKO30PP	DKO31PP	DKO32PP	DKO136PP		
24 mm cosmetic	3/8		45 blue		DKO154PP					
			75 blue		DKO143PP		DKO178PP			
26 mm	3/8		45 blue		DKO88PP	DKO36PP	DKO77PP			
			75 blue			DKO37PP	DKO38PP			
26 mm cosmetic	3/8		45 blue		DKO78PP	DKO39PP	DKO153PP			
			75 blue			DKO109PP	DKO108PP			
30 mm	3/8		45 blue			DKO47PP	DKO156PP			
			75 blue		DKO48PP	DKO49PP	DKO50PP	DKO90PP	DKO51PP	
30 mm cosmetic	3/8		75 blue				DKO158PP			
37 mm	3/8		75 blue				DKO61PP	DKO62PP		
			90 blue				DKO63PP			
39 mm	3/8		75 blue			DKO64PP	DKO65PP	DKO66PP	DKO67PP	
40 mm	3/8		45 blue				DKO126PP			
			75 blue					DKO105PP	DKO125PP	
			100 blue				DKO70PP		DKO144PP	
45 mm	3/8		75 blue				DKO71PP		DKO152PP	
			100 blue				DKO106PP			
77 mm	3/8		75 blue							DKO75PP
			3x60 blue							DKO95PP
77 mm	3/8		75 blue							DKO76PP

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



POLYPROPYLENE

surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
16 mm	1/2		45 blue		DKO180PP	DKO181PP				
16 mm	1/2		75 blue		DKO133PP					
19 mm	1/2		45 blue	DKO171PP	DKO170PP					
			75 blue		DKO163PP	DKO165PP	DKO149PP			
22 mm	1/2		45 blue			DKO23PP				
			75 blue			DKO24PP	DKO25PP			
22 mm	1/2		75 blue		DKO134PP					
24 mm	1/2		75 blue					DKO33PP		
			90 blue			DKO34PP				
26 mm	1/2		75 blue			DKO41PP	DKO42PP	DKO43PP		
			90 blue			DKO44PP	DKO45PP	DKO46PP		
			100 blue					DKO131PP	DKO132PP	
30 mm	1/2		75 blue				DKO150PP	DKO128PP	DKO129PP	
30 mm heavy	1/2		100 blue					DKO182PP	DKO121PP	
36 mm	1/2		75 blue				DKO151PP	DKO55PP	DKO56PP	DKO57PP
			90 blue					DKO58PP	DKO59PP	DKO60PP
40 mm	1/2		75 blue					DKO164PP		DKO84PP
			100 blue					DKO69PP		
48 mm	1/2		75 blue						DKO73PP	
			100 blue						DKO74PP	

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
51 mm	straight		90 blue				DG01PP			
55 mm	straight		100 blue				DG011PP			
60 mm	straight		75 blue			DG02PP	DG08PP			
60 mm	straight		75 blue				DG016PP			
70 mm	straight		90 blue					DG012PP		
75 mm	straight		100 blue				DG07PP			
77 mm	straight		75 blue				DG013PP			

notes



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



**cosmetic reverse cutting needle
with a precise blade (micro point)**

length	needles curvature		thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
13 mm	3/8		45 blue				ES10PP	ES08PP		
			75 blue					ES12PP		
16 mm	3/8		45 blue				ES11PP	ES03PP	ES04PP	
19 mm	3/8		45 blue					ES01PP	ES02PP	ES09PP
			75 blue						ES05PP	
26 mm	3/8		75 blue						ES06PP	ES07PP

notes



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



round needle with a cutting edge (sharpened)

length	needles curvature		thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
9 mm	3/8		75 blue		DS01PP	DS02PP				
9 mm	3/8		60 blue	DS37PP	DS33PP	DS46PP	DS57PP			
			75 blue		DS32PP	DS34PP				
9,3 mm	3/8		60 blue	DS107PP	DS94PP	DS88PP				
9,5 mm	3/8		60 blue	DS75PP						
10 mm	3/8		60 blue		DS100PP	DS99PP				
11 mm	3/8		60 blue			DS91PP				
11 mm	3/8		60 blue		DS59PP	DS78PP	DS47PP			
13 mm	3/8		75 blue			DS96PP				
13 mm	3/8		60 blue		DS08PP	DS38PP				
			75 blue			DS09PP	DS10PP			
17 mm	3/8		75 blue				DS97PP			
17 mm	3/8		90 blue				DS89PP	DS98PP	DS81PP	
20 mm	3/8		90 blue					DS121PP		
11 mm	1/2		60 blue			DS90PP	DS106PP			
13 mm	1/2		60 blue		DS66PP	DS53PP	DS54PP			
			75 blue				DS77PP	DS76PP		
16 mm	1/2		90 blue				DS44PP	DS45PP	DS68PP	

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



round needle with a cutting edge (sharpened)

length	needles curvature		thread length (cm) color	6/0	5/0	4/0	3/0	2/0	0	1
17 mm	1/2		75 blue			DS71PP				
			90 blue		DS48PP	DS13PP	DS14PP	DS15PP		
20 mm	1/2		90 blue			DS16PP	DS17PP			
22 mm	1/2		75 blue					DS118PP		
26 mm	1/2		75 blue			DS64PP	DS51PP	DS70PP	DS119PP	
			90 blue				DS21PP	DS36PP		
26 mm	1/2		90 blue		DS63PP	DS22PP	DS23PP	DS24PP		
			120 blue				DS112PP			
30 mm	1/2		90 blue				DS108PP	DS55PP		
			100 blue					DS69PP		
35 mm	1/2		75 blue						DS73PP	
			100 blue						DS25PP	
36 mm	1/2		90 blue						DS40PP	
36 mm	1/2		90 blue					DS26PP	DS79PP	

notes



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



round needle with a cutting edge (sharpened)

length	needles curvature		thread length (cm) color	6/0	5/0	4/0	3/0	2/0	0	1
37 mm	1/2		75 blue						DS27PP	
			90 blue							DS28PP
			100 blue						DS49PP	
40 mm	1/2		90 blue							DS74PP
40 mm	1/2		90 blue				DS87PP	DS58PP		
			100 blue						DS80PP	
45 mm	1/2		100 blue							DS31PP
48 mm	1/2		100 blue							DS116PP
30 mm	hook		100 blue						DS109PP	

notes



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



round cutting needle (sharpened)

length	needles curvature		thread length (cm) color	6/0	5/0	4/0	3/0	2/0	0	1
16 mm	3/8		75 blue						RD03PP	
22 mm	1/2		75 blue							RD04PP
25 mm	1/2		75 blue						RD05PP	RD06PP
25 mm	1/2		90 blue						RD01PP	
30 mm	1/2		75 blue							RD07PP RD08PP
36 mm	1/2		100 blue							RD02PP

notes



POLYPROPYLENE


surgical sutures, non-absorbable, monofilament



round needle with a trocar edge

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
22 mm	1/2		75 blue			TS03PP	TS04PP	TS06PP		
			90 blue				TS05PP			
22 mm	1/2		90 blue			TS17PP				
26 mm	1/2		75 blue				TS07PP	TS08PP		
			90 blue					TS09PP		
26 mm	1/2		90 blue			TS01PP	TS02PP			
30 mm	1/2		75 blue					TS10PP	TS11PP	
			90 blue						TS12PP	
36 mm	1/2		75 blue					TS13PP	TS14PP	TS16PP
			90 blue						TS15PP	
40 mm	1/2		75 blue					TS20PP	TS21PP	
48 mm	1/2		75 blue						TS18 PP	

notes



POLYPROPYLENE

surgical sutures, non-absorbable, monofilament



reverse cutting P needle blade

length	needles curvature		thread length (cm) color	7/0	6/0	5/0	4/0	3/0	2/0	0
13 mm	3/8		45 blue		DM02PP	DM05PP				
			75 blue	DM06PP						
16 mm	3/8		45 blue		DM03PP	DM01PP	DM04PP			
19 mm	3/8		75 blue					DM07PP		

notes

POLYPROPYLENE

surgical sutures, non-absorbable, monofilament

LIGATURES

threads length	color thread	7/0	6/0	5/0	4/0	3/0	2/0	0	1	2
2 x 100 cm	blue									P07PP
70 cm	blue				P01PP	P02PP	P03PP	P04PP		
75 cm	blue				P05PP					
150 cm	blue				P10PP	P11PP	P06PP	P08PP	P09PP	P12PP

notes