



PVDF

PVDF

Non-absorbable, synthetic surgical monofilament thread




- Polyvinylidene-fluoride
- Uncoated
- Thread may be blue or black
- Size 10/0 – 2 USP
- Features
 - Uniform structure with smooth surface
 - Very good binding properties
 - Tensile strength
- Indications
 - Soft tissues and skin
 - Traumatic surgery
 - Plastic surgery
 - Suturing and integration of blood vessels, also in cardiosurgery
 - Urology
 - Ophthalmology
- Period of validity 5 years

yavo

YAVO Sp. z o. o. 97-400 Bełchatów, ul. Bawełniana 17, Polska
tel.: +48 44 632 94 11, fax:+48 44 632 49 17,
e-mail: info@yavo.com.pl, www.yavo.com.pl



PVDF



PVDF (polyvinylidene-fluoride)


surgical sutures, non-absorbable, monofilament



spatula (lanced needle)

needles				thread length (cm) color	10/0	9/0	8/0	7/0	6/0	5/0	4/0
length	curvature	diameter									
6,19 mm	3/8	150 µm		30 blue	LZ10PV	LZ09PV					
6,4 mm	3/8	200 µm		30 blue		LZ07PV	LZ08PV				

notes



PVDF (polyvinylidene-fluoride)


surgical sutures, non-absorbable, monofilament



round needle

needles			thread length (cm) color							
length	curvature			9/0	8/0	7/0	6/0	5/0	4/0	3/0
6,4 mm	3/8		15 blue	RS02PV	RS03PV					
6,4 mm	3/8		45 blue	RS04PV	RS05PV					
8 mm	3/8		75 blue				RS64PV			
9,5 mm	3/8		75 blue			RS09PV				
11 mm	3/8		75 blue				RS20PV			
13 mm	3/8		75 blue				RS46PV	RS21PV		
16 mm	3/8		75 blue						RS34PV	
12 mm	1/2		75 blue				RS56PV	RS35PV		
16 mm	1/2		90 blue					RS47PV		
17 mm	1/2		75 blue					RS41PV	RS40PV	
22 mm	1/2		75 blue					RS62PV	RS42PV	RS17PV

notes



PVDF (polyvinylidene-fluoride)


surgical sutures, non-absorbable, monofilament



round needle

needles			thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
length	curvature									
22 mm	1/2		75 blue			RS16PV				
26 mm	1/2		75 blue		RS43PV	RS19PV	RS15PV	RS58PV		
26 mm	1/2		75 blue			RS25PV	RS26PV			
			90 blue			RS28PV	RS30PV			
30 mm	1/2		75 blue				RS14PV	RS54PV	RS60PV	
30 mm heavy	1/2		75 blue			RS61PV	RS13PV			
36 mm heavy	1/2		90 blue					RS18PV		
37 mm	1/2		75 blue				RS63PV	RS31PV	RS53PV	

notes



PVDF (polyvinylidene-fluoride)

surgical sutures, non-absorbable, monofilament



black needle

needles			thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
length	curvature									
12 mm	3/8		75 blue			RSC01PV				
17 mm	1/2		75 blue				RSC02PV			

notes



PVDF (polyvinylidene-fluoride)

surgical sutures, non-absorbable, monofilament



cutting needle

needles		thread length (cm) color								
length	curvature		6/0	5/0	4/0	3/0	2/0	0	1	
16 mm	3/8		45 blue		DK12PV					
			75 blue			DK06PV				
19 mm	3/8		45 blue			DK07PV				
20 mm	3/8		75 blue			DK03PV				
24 mm	3/8		45 blue			DK04PV				
			75 blue				DK08PV			
26 mm	3/8		75 blue				DK10PV	DK09PV	DK11PV	
30 mm	3/8		75 blue				DK02PV	DK01PV		
36 mm	3/8		75 blue				DK05PV			
60 mm	straight		90 blue					DG02PV		

notes



PVDF (polyvinylidene-fluoride)


surgical sutures, non-absorbable, monofilament



reverse cutting needle

length	needles curvature		thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
12 mm	3/8		45 blue				DKO30PV			
16 mm	3/8		45 blue				DKO18PV	DKO19PV	DKO36PV	
			75 blue				DKO20PV	DKO21PV	DKO32PV	DKO31PV
19 mm	3/8		45 blue				DKO09PV	DKO39PV	DKO28PV	DKO29PV
			75 blue						DKO44PV	DKO14PV
24 mm	3/8		45 blue						DKO43PV	DKO45PV
			75 blue							DKO46PV
			90 blue							DKO42PV
25 mm	3/8		45 blue							DKO13PV
			75 blue						DKO11PV	DKO15PV
30 mm	3/8		75 blue							DKO16PV
39 mm	3/8		75 blue							DKO41PV
8 mm	1/2		45 blue	DKO27PV						
60 mm	straight		75 blue							DKO1PV

notes



PVDF (polyvinylidene-fluoride)


surgical sutures, non-absorbable, monofilament



**cosmetic reverse cutting needle
with a precise blade (micro point)**

needles			thread length (cm) color	6/0	5/0	4/0	3/0	2/0	0	1
length	curvature									
13 mm	3/8		45 blue		ES09PV	ES01PV	ES02PV			
			75 blue		ES05PV					
16 mm	3/8		45 blue			ES04PV				
			75 blue			ES10PV				
18 mm	3/8		45 blue			ES07PV	ES06PV			
19 mm	3/8		45 blue		ES03PV					
			75 blue				ES11PV			
24 mm	3/8		75 blue			ES08PV		ES12PV		

notes



PVDF (polyvinylidene-fluoride)


surgical sutures, non-absorbable, monofilament



round needle with a cutting edge (sharpened)

needles		thread								
length	curvature		length (cm)	6/0	5/0	4/0	3/0	2/0	0	1
			color							
12,6 mm	3/8		60 blue	DS17PV						
			75 blue	DS25PV	DS14PV					
12,6 mm	3/8		75 blue	DS06PV	DS07PV					
16 mm	1/2		75 blue		DS24PV					
17 mm	1/2		75 blue		DS10PV		DS09PV			
17 mm	1/2		90 blue			DS15PV				
20 mm	1/2		75 blue					DS08PV		
25 mm	1/2		90 blue					DS20PV		
26 mm	1/2		75 blue				DS22PV	DS21PV		
30 mm	1/2		90 blue					DS13PV		
35 mm	1/2		75 blue						DS12PV	
35 mm	1/2		90 blue				DS16PV			

notes



PVDF (polyvinylidene-fluoride)

surgical sutures, non-absorbable, monofilament



reverse cutting P needle blade

length	needles		thread	6/0	5/0	4/0	3/0	2/0	0	1
	curvature		length (cm)							
	color									
13 mm	3/8		45 blue		DM04PV					
16 mm	3/8		45 blue			DM01PV				
18 mm	3/8		45 blue			DM03PV	DM02PV			

notes