



**PGA RAPID**


# PGA RAPID

Absorbable, synthetic, braided surgical thread

Uniform structure, strong tightening of knots, minimum tissue response, exceptional structural strength.



- 100% polyglycolid acid
- Polycaprolactone and calcium stearate coated
- Complete absorption after 42 days
- Thread undyed
- Size 6/0 – 2 USP
- Breaking strenght retention




- Indications
  - Soft tissues
  - Closing the dermal wounds, especially in pediatric surgery
  - Joining mucous membrane in mouth
  - General and plastic surgery
  - Gynecology
  - Ophthalmic surgery
- Period of validity 3 years

**yavo**

YAVO Sp. z o. o. 97-400 Bełchatów, ul. Bawełniana 17, Polska  
tel.: +48 44 632 94 11, fax:+48 44 632 49 17,  
e-mail: info@yavo.com.pl, www.yavo.com.pl



PGA RAPID
















## PGA RAPID (polyglycolid acid)


surgical sutures, absorbable about 42 days,  
multifilament, coated



**round needle**

needles		thread length (cm) color	6/0	5/0	4/0	3/0	2/0	0	1
length	curvature								
16 mm	3/8	 75 undyed		RS59PR					
22 mm	3/8	 90 undyed				RS109PR			
26 mm	3/8	 75 undyed			RS132PR				
							RS110PR		
30 mm	3/8	 75 undyed					RS60PR		
13 mm	1/2	 45 undyed		RS65PR					
				RS01PR					
16 mm	1/2	 75 undyed	RS55PR	RS02PR	RS03PR	RS04PR			
17 mm	1/2	 45 undyed			RS98PR	RS123PR			
				RS78PR	RS69PR	RS90PR			
20 mm	1/2	 75 undyed			RS87PR	RS05PR	RS86PR		
22 mm	1/2	 75 undyed		RS06PR	RS07PR	RS08PR	RS09PR		
						RS63PR	RS76PR		
26 mm	1/2	 45 undyed			RS99PR	RS83PR			
					RS10PR	RS11PR	RS12PR	RS13PR	
						RS111PR	RS62PR		
26 mm heavy	1/2	 70 undyed					RS93PR	RS94PR	
								RS50PR	
26 mm thin	1/2	 75 undyed					RS52PR		
30 mm	1/2	 75 undyed			RS115PR	RS113PR	RS15PR	RS16PR	RS17PR
							RS64PR	RS106PR	

notes



## PGA RAPID (polyglycolid acid)

surgical sutures, absorbable about 42 days,  
multifilament, coated



round needle

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
30 mm heavy	1/2		90 undyed				RS18PR	RS19PR	RS20PR	
35 mm	1/2		75 undyed				RS82PR	RS117PR	RS118PR	RS119PR
			90 undyed				RS57PR	RS45PR	RS46PR	
37 mm	1/2		75 undyed			RS22PR	RS23PR	RS24PR	RS25PR	RS58PR
			90 undyed			RS26PR	RS27PR	RS28PR	RS29PR	RS54PR
37 mm heavy	1/2		75 undyed			RS127PR	RS30PR	RS31PR	RS32PR	
			90 undyed				RS33PR	RS34PR	RS35PR	RS48PR
40 mm	1/2		75 undyed				RS105PR	RS114PR	RS66PR	RS124PR
			90 undyed				RS56PR	RS36PR	RS37PR	
40 mm heavy	1/2		75 undyed						RS89PR	
			90 undyed					RS38PR	RS39PR	RS53PR
45 mm	1/2		75 undyed					RS67PR		RS125PR
			90 undyed				RS61PR	RS40PR	RS41PR	RS70PR
48 mm	1/2		75 undyed						RS100PR	RS120PR
			90 undyed				RS49PR	RS42PR	RS43PR	RS44PR
48 mm heavy	1/2		90 undyed				RS79PR	RS91PR	RS92PR	RS47PR
50 mm	1/2		90 undyed							RS126PR
76 mm	1/2		75 undyed				RS97PR			
27 mm	5/8		75 undyed				RS77PR			

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



## PGA RAPID (polyglycolid acid)

surgical sutures, absorbable about 42 days,  
multifilament, coated



### cutting needle

needles			thread length (cm) color	6/0	5/0	4/0	3/0	2/0	0	1
length	curvature									
12 mm	3/8		45 undyed		DK20PR					
16 mm	3/8		75 undyed		DK05PR	DK06PR	DK07PR			
16 mm cosmetic	3/8		45 undyed			DK03PR	DK04PR			
			75 undyed			DK01PR	DK02PR			
19 mm	3/8		75 undyed			DK08PR	DK09PR	DK10PR		
19 mm cosmetic	3/8		75 undyed				DK35PR			
24 mm	3/8		75 undyed				DK17PR	DK25PR		
26 mm	3/8		75 undyed					DK19PR		
26 mm cosmetic	3/8		75 undyed					DK26PR		
30 mm	3/8		75 undyed					DK18PR		
35 mm	3/8		75 undyed					DK21PR		
36 mm	1/2		75 undyed					DK36PR	DK37PR	
51 mm	straight		75 undyed					DK30PR		
60 mm	straight		75 undyed				DK27PR	DK33PR	DK34PR	



## PGA RAPID (polyglycolilid acid)

surgical sutures, absorbable about 42 days,  
multifilament, coated



## reverse cutting needle

needles			thread length (cm) color								
length	curvature			6/0	5/0	4/0	3/0	2/0	0	1	
12 mm	3/8		45 undyed	DKO01PR	DKO02PR						
12 mm cosmetic	3/8		45 undyed	DKO03PR	DKO04PR						
13 mm	3/8		45 undyed			DKO77PR					
			75 undyed		DKO78PR	DKO95PR					
13 mm cosmetic	3/8		45 undyed		DKO66PR	DKO67PR					
16 mm	3/8		45 undyed		DKO97PR	DKO73PR					
			75 undyed		DKO05PR	DKO06PR	DKO07PR				
16 mm cosmetic	3/8		45 undyed	DKO107PR	DKO106PR	DKO103PR					
			75 undyed		DKO08PR	DKO09PR	DKO10PR				
19 mm	3/8		45 undyed		DKO87PR	DKO18PR	DKO19PR	DKO20PR			
			75 undyed		DKO92PR	DKO15PR	DKO16PR	DKO17PR			
19 mm cosmetic	3/8		45 undyed			DKO102PR	DKO101PR				
			75 undyed			DKO13PR	DKO14PR				
24 mm	3/8		75 undyed		DKO70PR	DKO21PR	DKO22PR	DKO23PR			
24 mm cosmetic	3/8		45 undyed			DKO27PR	DKO28PR	DKO29PR			
			75 undyed			DKO24PR	DKO25PR	DKO26PR			
26 mm	3/8		45 undyed				DKO36PR	DKO37PR			
			75 undyed				DKO34PR	DKO35PR	DKO69PR		
26 mm cosmetic	3/8		45 undyed				DKO32PR	DKO33PR			
			75 undyed				DKO30PR	DKO31PR			
30 mm	3/8		75 undyed					DKO109PR	DKO93PR		
			90 undyed				DKO68PR	DKO38PR	DKO39PR		

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



## PGA RAPID (polyglycolid acid)


surgical sutures, absorbable about 42 days,  
multifilament, coated



## reverse cutting needle

length	needles curvature	thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
35 mm	3/8	75 undyed			DKO40PR	DKO41PR			
35 mm cosmetic	3/8	75 undyed			DKO42PR	DKO43PR			
		90 undyed			DKO44PR	DKO45PR	DKO46PR		
35 mm heavy	3/8	75 undyed				DKO64PR			
16 mm	1/2	75 undyed	DKO122PR	DKO100PR					
22 mm	1/2	45 undyed		DKO117PR					
24 mm	1/2	45 undyed		DKO89PR	DKO80PR				
		75 undyed			DKO86PR				
37 mm	1/2	75 undyed						DKO121PR	
		90 undyed				DKO54PR	DKO55PR	DKO56PR	
37 mm heavy	1/2	75 undyed				DKO75PR			
		90 undyed				DKO81PR			
40 mm	1/2	75 undyed					DKO96PR	DKO79PR	
		90 undyed						DKO57PR	
40 mm heavy	1/2	90 undyed						DKO65PR	
45 mm	1/2	90 undyed						DKO58PR	
48 mm	1/2	90 undyed					DKO63PR	DKO59PR	
51 mm	straight	70 undyed		DKO83PR	DKO84PR				
60 mm	straight	75 undyed			DKO94PR				

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



## PGA RAPID (polyglycolild acid)


surgical sutures, absorbable about 42 days,  
multifilament, coated



**cosmetic reverse cutting needle  
with a precise blade (micro point)**

needles			thread length (cm) color								
length	curvature			6/0	5/0	4/0	3/0	2/0	0	1	
16 mm	3/8		45 undyed		ES03PR						
19 mm	3/8		45 undyed		ES04PR	ES05PR	ES10PR				
			75 undyed			ES02PR					
24 mm	3/8		45 undyed			ES06PR					
			75 undyed			ES07PR					
26 mm	3/8		45 undyed			ES08PR	ES09PR				
			75 undyed				ES01PR				

notes



## PGA RAPID (polyglycolid acid)

surgical sutures, absorbable about 42 days,  
multifilament, coated




**round needle with a cutting  
edge (sharpened)**

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
15 mm	3/8		75 undyed	DS32PR						
17 mm	1/2		90 undyed		DS36PR					
20 mm	1/2		75 undyed		DS39PR					
22 mm	1/2		75 undyed			DS40PR				
30 mm	1/2		75 undyed				DS45PR	DS44PR		
			90 undyed				DS21PR			
35 mm	1/2		75 undyed				DS07PR	DS08PR	DS09PR	
			90 undyed				DS10PR	DS11PR	DS12PR	
37 mm	1/2		90 undyed				DS13PR	DS14PR	DS15PR	DS22PR
40 mm	1/2		75 undyed						DS51PR	
			90 undyed					DS35PR	DS16PR	
40 mm heavy	1/2		90 undyed						DS47PR	
45 mm	1/2		90 undyed					DS30PR	DS17PR	DS20PR
			120 undyed						DS42PR	
48 mm	1/2		90 undyed				DS29PR	DS18PR	DS19PR	
48 mm heavy	1/2		90 undyed					DS56PR		
17 mm	5/8		75 undyed		DS57PR					

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



## PGA RAPID (polyglycolid acid)


surgical sutures, absorbable about 42 days,  
multifilament, coated



**round needle with a trocar edge**

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
22 mm	1/2		75 undyed			TS01PR	TS02PR			
			90 undyed			TS03PR	TS04PR			
25 mm	1/2		75 undyed			TS05PR	TS06PR	TS07PR	TS08PR	
			90 undyed				TS09PR	TS10PR	TS11PR	
30 mm	1/2		75 undyed					TS12PR	TS13PR	
			90 undyed					TS14PR	TS15PR	
36 mm	1/2		75 undyed					TS16PR	TS17PR	TS18PR
			90 undyed				TS28PR	TS19PR	TS20PR	TS21PR
40 mm	1/2		75 undyed					TS22PR	TS23PR	TS24PR
			90 undyed					TS25PR	TS26PR	TS27PR
40 mm heavy	1/2		90 undyed					TS32PR	TS31PR	
45 mm heavy	1/2		90 undyed					TS30PR		

notatki




## PGA RAPID (polyglycolid acid)

surgical sutures, absorbable about 42 days,  
multifilament, coated



### round cutting needle (sharpened)

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
37 mm	1/2		90 undyed					RD03PR	RD02PR	
40 mm	1/2		90 undyed						RD01PR	



### round separating needle

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
17 mm	1/2		75 undyed		RR03PR					
22 mm	1/2		75 undyed			RR01PR				
26 mm	1/2		75 undyed		RR02PR					



### reverse cutting P needle blade

length	needles curvature		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
16 mm	3/8		75 undyed	DM06PR	DM01PR	DM02PR				
19 mm	3/8		45 undyed	DM05PR	DM04PR	DM03PR				
			75 undyed	DM09PR	DM08PR	DM07PR				

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.

## PGA RAPID (polyglycolild acid)

surgical sutures, absorbable about 42 days,  
multifilament, coated

### LIGATURES

threads length	color thread	7/0	6/0	5/0	4/0	3/0	2/0	0	1	2
2 x 75 cm	undyed					P16PR	P17PR	P18PR	P19PR	P20PR
3 x 45 cm	undyed				P07PR	P08PR	P09PR	P10PR	P11PR	P12PR
5 x 75 cm	undyed				P27PR	P28PR	P29PR	P30PR	P31PR	P32PR
150 cm	undyed				P41PR	P42PR	P43PR	P44PR	P45PR	P46PR
250 cm	undyed					P61PR	P47PR	P48PR	P49PR	P50PR
250 cm in reel	undyed						Pr47PR			

notes