


PGLA LACTIC


PGLA LACTIC

Absorbable, synthetic, braided surgical thread

Uniform structure, strong tightening of knots, minimum tissue response, exceptional structural strength



- Poly (glycolide-co- L- lactide)(90/10)
- Mixture of poly(glycolide-co-L-Lactide) (30/70) (50%) and calcium stearate (50%)
- Complete absorption after 60-70 days
- Thread may be natural (beige) or violet
- Size 8/0 – 2 USP
- Breaking strenght retention



- Indications
 - Soft tissues (subcutaneous tissues, mucous membrane) and/or tying off tissues
 - Pediatric, general surgery, plastic surgery
 - Gynecology/obstetrics
 - Ophthalmic surgery
 - Urology
 - Orthopedics
 - Dental surgery
- Period of validity 3 years

yavo

YAVO Sp. z o. o. 97-400 Bełchatów, ul. Bawełniana 17, Polska
tel.: +48 44 632 94 11, fax:+48 44 632 49 17,
e-mail: info@yavo.com.pl, www.yavo.com.pl



PGLA LACTIC



PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated



spatula (lanced needle)

needles		thread length (cm) color		8/0	7/0	6/0	5/0	4/0	3/0	2/0
length	curvature									
6 mm	3/8	230 µm		30 violet	LZ51PL	LZ43PL				
				45 violet	LZ01PL	LZ02PL	LZ58PL			
6 mm	3/8	230 µm		30 violet	LZ27PL	LZ28PL				
				45 violet	LZ03PL	LZ04PL	LZ05PL			
6 mm micropoint	3/8	220 µm		45 violet	LZ59PL					
6 mm micropoint	3/8	230 µm		30 violet	LZ48PL					
				45 violet		LZ47PL				
6,19 mm	3/8	230 µm		45 violet	LZ54PL					
6,19 mm	3/8	150 µm		45 violet	LZ77PL					
6,4 mm	3/8	230 µm		20 violet		LZ53PL				
				45 violet	LZ17PL	LZ18PL				
6,5 mm micropoint	3/8	220 µm		30 violet	LZ78PL	LZ68PL				
				45 violet	LZ69PL					
6,5 mm	3/8	230 µm		30 violet	LZ20PL	LZ16PL				
				45 violet	LZ42PL					
6,5 mm micropoint	3/8	230 µm		30 violet	LZ52PL					
				45 violet	LZ61PL					
8 mm	3/8	230 µm		45 violet			LZ25PL			
10 mm	3/8	230 µm		45 violet			LZ13PL			
6 mm	1/2			30 violet		LZ75PL				
6,5 mm	1/2	250 µm		30 violet	LZ34PL					
				45 violet	LZ39PL					
8 mm	1/2			45 undyed			LZ70PL			
				45 violet		LZ80PL				

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



spatula (lanced needle)

needles			thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
length	curvature	diameter								
8 mm	1/2	300 µm		45 undyed			LZ63PL			
				45 violet	LZ09PL	LZ10PL	LZ11PL	LZ12PL		
8 mm micropoint	1/2	220 µm		45 violet		LZ60PL	LZ71PL			
8 mm micropoint	1/2	300 µm		45 violet	LZ56PL					
8 mm	1/4	350 µm		20 violet			LZ21PL			
				45 violet			LZ06PL	LZ07PL	LZ08PL	
				45 undyed			LZ31PL	LZ26PL		
				6x45 violet			LZ41PL			
8 mm micropoint	1/4	330 µm		45 violet			LZ73PL	LZ62PL		
				45 undyed				LZ72PL		

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round needle

needles		thread length (cm) color								
length	curvature		8/0	7/0	6/0	5/0	4/0	3/0	2/0	
6 mm	3/8		13 violet	RS281PL						
			30 violet	RS385PL						
8 mm	3/8		30 violet	RS350PL						
			45 violet	RS01PL	RS02PL	RS03PL				
10 mm	3/8		45 violet		RS493PL					
13 mm	3/8		45 violet			RS09PL	RS10PL	RS11PL		
			60 violet					RS221PL		
			75 violet			RS06PL	RS07PL	RS08PL		
17 mm	3/8		45 violet					RS268PL		
			75 violet				RS322PL	RS199PL		
22 mm	3/8		45 violet				RS411PL			
			75 violet					RS321PL	RS166PL	RS159PL
26 mm	3/8		75 violet						RS168PL	RS49PL
			75 undyed							RS224PL

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round needle

needles		thread length (cm) color								
length	curvature		6/0	5/0	4/0	3/0	2/0	0	1	
32 mm	3/8		75 violet				RS215PL	RS360PL	RS286PL	RS469PL
13 mm	1/2		45 violet	RS15PL	RS16PL	RS17PL				
			45 undyed	RS328PL	RS330PL					
			60 violet			RS246PL				
			75 violet	RS12PL	RS13PL	RS14PL				
			75 undyed		RS391PL	RS426PL				
16 mm	1/2		45 violet		RS21PL	RS22PL				
			75 violet		RS19PL	RS20PL	RS153PL	RS174PL		
			75 undyed		RS274PL					
16 mm thin	1/2		75 violet				RS461PL			
17 mm	1/2		18 violet				RS476PL			
			45 undyed			RS512PL				
			75 violet	RS173PL	RS162PL	RS163PL	RS181PL	RS203PL		
			75 undyed		RS511PL	RS278PL				
17 mm thin	1/2		75 violet				RS450PL			
20 mm	1/2		45 violet			RS27PL	RS28PL	RS272PL		
			75 violet		RS225PL	RS25PL	RS26PL	RS190PL		
			75 undyed			RS275PL				
20 mm thin	1/2		75 violet			RS342PL	RS340PL			
22 mm thin	1/2		75 violet			RS193PL	RS367PL			
22 mm	1/2		90 violet			RS32PL	RS33PL	RS34PL		

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.





PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round needle

needles		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2	
length	curvature									
22 mm	1/2		45 violet		RS245PL	RS241PL	RS235PL			
		75 violet		RS29PL	RS30PL	RS31PL	RS197PL		RS409PL	
		75 undyed		RS503PL	RS390PL	RS502PL				
26 mm	1/2		45 violet			RS243PL	RS238PL	RS229PL		
		75 violet	RS439PL	RS37PL	RS38PL	RS39PL	RS40PL	RS41PL	RS42PL	
		75 undyed		RS331PL	RS276PL	RS277PL				
		90 violet			RS218PL	RS160PL	RS171PL		RS497PL	
26 mm thin	1/2		75 violet		RS352PL	RS155PL	RS156PL	RS372PL		
26 mm heavy	1/2		75 violet			RS332PL	RS161PL	RS169PL		RS507PL
			90 violet					RS380PL		RS516PL
26 mm	1/2		90 violet			RS44PL	RS45PL	RS46PL		
30 mm	1/2		75 violet		RS373PL	RS55PL	RS56PL	RS57PL	RS58PL	RS265PL
			75 undyed				RS456PL			
			90 violet			RS437PL	RS269PL	RS359PL	RS259PL	RS320PL
			100 violet				RS412PL			
30 mm thin	1/2		75 violet			RS51PL	RS52PL	RS53PL	RS54PL	
30 mm heavy	1/2		75 violet			RS59PL	RS60PL	RS61PL	RS62PL	RS227PL
			90 violet						RS418PL	
30 mm	1/2		75 violet				RS64PL			
			100 violet					RS250PL		

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round needle

needles		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2	
length	curvature									
35 mm	1/2		75 violet				RS67PL	RS68PL	RS214PL	RS206PL
			90 violet				RS65PL	RS66PL	RS220PL	RS515PL
35 mm heavy	1/2		75 violet				RS434PL	RS362PL	RS363PL	
			90 violet				RS366PL	RS436PL		RS200PL
37 mm	1/2		45 violet					RS228PL	RS230PL	
			75 violet		RS432PL	RS71PL	RS72PL	RS73PL	RS74PL	RS75PL
			75 undyed			RS346PL	RS347PL			
			90 violet			RS76PL	RS77PL	RS78PL	RS79PL	RS80PL
			90 undyed						RS446PL	
			120 violet				RS195PL			
37 mm heavy	1/2		75 violet			RS81PL	RS82PL	RS83PL	RS84PL	RS85PL
			90 violet				RS451PL	RS337PL	RS204PL	RS338PL
37 mm thin	1/2		75 violet			RS341PL	RS247PL	RS361PL	RS479PL	
40 mm	1/2		75 violet			RS211PL	RS93PL	RS94PL	RS95PL	RS96PL
			90 violet			RS287PL	RS89PL	RS90PL	RS91PL	RS92PL
40 mm heavy	1/2		75 violet				RS97PL	RS98PL	RS99PL	RS100PL
			90 violet				RS500PL	RS217PL	RS150PL	RS290PL
40 mm	1/2		150 violet					RS105PL	RS106PL	
40 mm heavy	1/2		150 violet					RS495PL		

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated




round needle

needles		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
length	curvature								
45 mm	1/2	75 violet				RS107PL	RS108PL	RS109PL	RS110PL
		90 violet				RS111PL	RS112PL	RS113PL	RS114PL
45 mm heavy	1/2	75 violet					RS115PL	RS116PL	RS117PL
		90 violet							RS427PL
48 mm	1/2	75 violet			RS212PL	RS118PL	RS119PL	RS120PL	RS121PL
		90 violet			RS242PL	RS122PL	RS123PL	RS124PL	RS125PL
48 mm heavy	1/2	75 violet			RS378PL	RS379PL	RS126PL	RS127PL	RS128PL
		90 violet					RS381PL	RS263PL	RS283PL
48 mm	1/2	150 violet loop						RS498PL	
50 mm	1/2	75 violet					RS135PL	RS136PL	RS137PL
		90 violet					RS132PL	RS133PL	RS134PL
50 mm heavy	1/2	90 violet						RS387PL	RS388PL
60 mm	1/2	75 violet				RS138PL	RS139PL	RS140PL	RS405PL
60 mm	1/2	150 violet loop						RS440PL	
60 mm heavy	1/2	90 violet							RS431PL

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round needle

needles		thread		5/0	4/0	3/0	2/0	0	1	2
length	curvature	length (cm)	color							
65 mm	1/2	75	violet				RS158PL		RS231PL	
		75	undyed				RS357PL			
		90	violet				RS141PL	RS142PL	RS143PL	RS144PL
		150	violet						RS170PL	RS154PL
65 mm heavy	1/2	75	violet				RS288PL		RS370PL	RS301PL
		90	violet						RS291PL	RS327PL
		150	violet						RS485PL	
65 mm	1/2	150	violet						RS483PL	
		240	violet							RS267PL
76 mm	1/2	150	violet							RS480PL
76 mm	1/2	75	violet			RS244PL	RS237PL	RS398PL	RS399PL	RS406PL
		75	undyed				RS445PL			
		90	violet						RS371PL	RS441PL
76 mm heavy	1/2	75	violet				RS463PL		RS482PL	
		90	violet							RS393PL
		100	violet						RS410PL	
26 mm	5/8	75	violet			RS423PL	RS306PL	RS305PL		
27 mm	5/8	75	violet			RS280PL	RS50PL	RS172PL	RS165PL	

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated




round needle

needles			thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
length	curvature									
27 mm heavy	5/8		75 violet						RS300PL	
36 mm	5/8		75 violet				RS70PL	RS317PL	RS448PL	
37 mm	5/8		75 violet				RS202PL	RS266PL		
40 mm	5/8		75 violet					RS180PL	RS468PL	
			90 violet					RS514PL		
22 mm	hook		75 violet		RS339PL	RS35PL				
30 mm	hook		75 violet				RS316PL	RS315PL	RS63PL	
31 mm	hook		75 violet						RS304PL	
32 mm	hook		75 violet						RS402PL	
35 mm	hook		75 violet				RS324PL	RS69PL	RS449PL	
37 mm	hook		70 violet					RS443PL		
			75 violet						RS408PL	

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round needle

needles			thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
length	curvature									
22 mm	ski		20 violet				RS492PL			
23 mm	ski		18 violet		RS491PL					
			20 violet			RS403PL				
			100 violet			RS435PL				
			120 violet			RS478PL				
23 mm	ski	 loop	20 violet			RS36PL	RS513PL			
26 mm	ski		120 violet			RS505PL	RS506PL			
19 mm	straight		20 violet			RS23PL				

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round blunt needle

needles		thread length (cm) color							
length	curvature		5/0	4/0	3/0	2/0	0	1	2
63 mm	3/8	75 violet						RST02PL	
		100 violet						RST03PL	
65 mm	3/8	100 violet						RST04PL	
30 mm	1/2	75 violet			RST05PL	RST06PL			
35 mm	1/2	75 violet				RST08PL			
35 mm heavy	1/2	90 violet					RST09PL	RST10PL	
37 mm	1/2	75 violet				RST12PL	RST13PL	RST14PL	
		75 undyed				RST15PL			
		90 violet					RST16PL	RST17PL	RST18PL
		90 undyed			RST19PL	RST20PL			
40 mm	1/2	75 violet					RST23PL	RST24PL	
		90 violet				RST25PL	RST26PL	RST27PL	RST28PL
48 mm	1/2	90 violet					RST32PL	RST33PL	RST34PL
48 mm heavy	1/2	90 violet							RST35PL
65 mm	1/2	75 violet							RST37PL
		100 violet						RST38PL	RST39PL
		150 violet						RST41PL	RST42PL

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated



round black needle

needles			thread length (cm) color							
length	curvature			6/0	5/0	4/0	3/0	2/0	0	1
13 mm	3/8		75 violet	RSC01PL						
23 mm	1/2		75 violet					RSC03PL		
26 mm	1/2		75 violet				RSC04PL			

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated



cutting needle

needles			thread length (cm) color	6/0	5/0	4/0	3/0	2/0	0	1
length	curvature									
12 mm	3/8		45 violet	DK109PL						
			45 undyed	DK146PL						
12 mm cosmetic	3/8		45 violet		DK65PL					
			45 undyed		DK137PL					
16 mm	3/8		45 violet	DK89PL	DK94PL	DK67PL	DK68PL			
			75 violet		DK10PL	DK11PL	DK12PL			
			75 undyed			DK147PL	DK129PL			
16 mm cosmetic	3/8		45 violet	DK135PL	DK07PL	DK08PL	DK09PL			
			45 undyed		DK138PL	DK139PL				
19 mm	3/8		45 violet		DK23PL	DK24PL	DK25PL			
			45 undyed			DK110PL	DK83PL			
			75 violet		DK20PL	DK21PL	DK22PL			
			75 undyed		DK149PL	DK125PL	DK124PL			
19 mm cosmetic	3/8		45 violet		DK17PL	DK18PL	DK19PL			
			75 violet		DK14PL	DK15PL	DK16PL			
24 mm	3/8		45 violet			DK32PL	DK33PL	DK34PL		
			45 undyed				DK112PL			
			75 violet			DK29PL	DK30PL	DK31PL		
			75 undyed			DK95PL	DK96PL	DK97PL		
26 mm	3/8		75 violet			DK41PL	DK42PL	DK43PL	DK44PL	
26 mm cosmetic	3/8		45 violet			DK38PL	DK39PL	DK40PL		
			75 undyed			DK140PL				

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.




PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated



cutting needle

needles		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
length	curvature								
30 mm	3/8	 75 violet			DK70PL	DK103PL	DK98PL		
36 mm	3/8	 75 violet			DK81PL	DK79PL			
36 mm cosmetic	3/8	 75 violet			DK78PL	DK101PL			
		75 undyed			DK104PL	DK145PL			
40 mm	3/8	 75 violet						DK128PL	
22 mm	1/2	 45 violet			DK87PL				
		75 violet			DK148PL	DK144PL			DK108PL
24 mm	1/2	 45 undyed			DK111PL				
		75 undyed		DK93PL	DK90PL				
26 mm	1/2	 45 violet			DK132PL				
		75 violet					DK69PL		DK102PL
37 mm	1/2	 75 violet			DK53PL	DK54PL	DK55PL	DK56PL	DK57PL
		90 violet				DK88PL	DK85PL		
40 mm	1/2	 75 violet				DK58PL	DK59PL	DK60PL	DK61PL
		90 violet						DK76PL	DK66PL
13 mm	straight	 45 violet		DG01PL					
45 mm	straight	 45 violet		DG02PL					
55 mm	straight	 75 violet				DG03PL			
60 mm	straight	 75 violet			DG04PL	DG05PL	DG08PL		
		75 undyed			DG07PL	DG06PL			

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated



reverse cutting needle

needles			thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
length	curvature									
9 mm	3/8		45 violet			DKO145PL				
12 mm	3/8		45 violet			DKO03PL	DKO04PL	DKO05PL		
			45 undyed				DKO190PL			
12 mm cosmetic	3/8		45 violet				DKO102PL			
			45 undyed			DKO186PL				
13 mm	3/8		45 violet			DKO119PL	DKO81PL			
16 mm	3/8		45 violet			DKO125PL	DKO06PL	DKO07PL	DKO08PL	
			75 violet				DKO09PL	DKO10PL	DKO11PL	
			75 undyed				DKO210PL		DKO176PL	
16 mm cosmetic	3/8		45 violet				DKO12PL	DKO13PL	DKO14PL	
			45 undyed			DKO189PL	DKO161PL	DKO133PL		
19 mm	3/8		45 violet				DKO18PL	DKO19PL	DKO20PL	
			45 undyed				DKO150PL	DKO122PL	DKO198PL	
			75 violet				DKO67PL	DKO21PL	DKO22PL	DKO147PL
19 mm cosmetic	3/8		45 violet				DKO15PL	DKO16PL	DKO17PL	
			45 undyed				DKO201PL	DKO128PL	DKO121PL	
			75 undyed					DKO168PL		
24 mm	3/8		45 violet					DKO85PL	DKO204PL	DKO75PL
			75 violet					DKO26PL	DKO27PL	DKO28PL

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated



reverse cutting needle

needles		thread length (cm) color							
length	curvature		5/0	4/0	3/0	2/0	0	1	2
26 mm	3/8	45 violet			DKO167PL	DKO34PL			
		75 violet			DKO158PL	DKO33PL	DKO66PL	DKO179PL	
		75 undyed			DKO129PL	DKO91PL	DKO101PL		
26 mm cosmetic	3/8	45 violet				DKO31PL	DKO32PL		
		75 violet				DKO29PL	DKO30PL		
		75 undyed					DKO188PL		
30 mm	3/8	75 violet				DKO63PL	DKO65PL	DKO68PL	
35 mm	3/8	75 violet				DKO39PL	DKO40PL	DKO41PL	DKO72PL
37 mm	3/8	75 violet					DKO149PL		
37 mm cosmetic	3/8	75 undyed				DKO113PL			
40 mm	3/8	75 violet						DKO73PL	
		90 violet							DKO181PL
45 mm	3/8	75 violet					DKO89PL		
48 mm	3/8	90 violet						DKO208PL	
60 mm	3/8	90 violet						DKO192PL	

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days, multifilament, coated



reverse cutting needle

length	needles		thread length (cm) color	8/0	7/0	6/0	5/0	4/0	3/0	2/0
	curvature									
8 mm	1/2		30 violet	DKO64PL						
16 mm	1/2		45 undyed					DKO196PL		
19 mm	1/2		75 violet				DKO177PL	DKO216PL	DKO76PL	
22 mm	1/2		45 violet						DKO165PL	

notes





PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



reverse cutting needle

needles		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2	
length	curvature									
22 mm heavy	1/2		75 violet							DKO99PL
24 mm	1/2		75 violet			DKO77PL	DKO104PL			
26 mm	1/2		75 violet			DKO182PL	DKO70PL	DKO88PL	DKO69PL	DKO35PL
26 mm heavy	1/2		75 violet				DKO78PL	DKO185PL	DKO152PL	DKO184PL
37 mm	1/2		75 violet			DKO42PL	DKO43PL	DKO44PL	DKO45PL	DKO62PL
			90 violet					DKO97PL		DKO94PL
37 mm heavy	1/2		75 violet				DKO79PL		DKO112PL	DKO162PL
			90 violet						DKO195PL	
40 mm	1/2		75 violet			DKO134PL	DKO49PL	DKO50PL	DKO51PL	DKO172PL
			90 violet				DKO46PL	DKO47PL	DKO48PL	DKO71PL
40 mm heavy	1/2		90 violet				DKO193PL			
45 mm	1/2		75 violet				DKO52PL		DKO87PL	
48 mm	1/2		75 violet				DKO154PL	DKO153PL	DKO156PL	DKO123PL
			90 violet					DKO207PL	DKO105PL	DKO92PL
48 mm heavy	1/2		90 violet					DKO138PL	DKO157PL	DKO215PL
50 mm	1/2		75 violet					DKO56PL	DKO57PL	DKO58PL
			90 violet					DKO53PL	DKO54PL	DKO55PL
60 mm	straight		75 violet				DG09PL			

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



**cosmetic reverse cutting needle
with a precise blade (micro point)**

needles		thread length (cm) color							
length	curvature		6/0	5/0	4/0	3/0	2/0	0	1
13 mm	3/8		45 violet	ES19PL	ES13PL	ES12PL			
			45 undyed	ES20PL	ES11PL				
16 mm	3/8		45 violet			ES17PL	ES18PL		
			45 undyed		ES01PL				
19 mm	3/8		45 violet		ES03PL	ES04PL			
			45 undyed			ES14PL	ES10PL		
			75 violet			ES05PL			
24 mm	3/8		75 violet			ES06PL	ES07PL		
			75 undyed				ES15PL		
26 mm	3/8		45 violet			ES21PL	ES22PL	ES23PL	
			45 undyed			ES08PL	ES02PL		
			75 violet				ES09PL	ES16PL	
30 mm	3/8		75 violet					ES24PL	

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated




**round needle with a cutting
edge (sharpened)**

needles			thread length (cm) color							
length	curvature			5/0	4/0	3/0	2/0	0	1	2
13 mm	3/8		45 violet	DS04PL	DS05PL					
16 mm	3/8		75 violet	DS06PL	DS07PL	DS08PL				
17 mm	3/8		45 violet		DS147PL					
			75 violet		DS54PL	DS51PL				
22 mm	3/8		75 violet		DS99PL		DS95PL			
24 mm	3/8		75 violet		DS13PL	DS14PL	DS15PL			
17 mm	1/2		75 violet	DS114PL	DS131PL	DS158PL				
20 mm	1/2		75 violet	DS143PL	DS09PL	DS10PL	DS11PL			
22 mm	1/2		75 violet		DS112PL	DS105PL	DS52PL	DS125PL		
26 mm	1/2		75 violet		DS111PL	DS16PL	DS17PL	DS85PL	DS61PL	
			90 violet				DS141PL		DS47PL	
26 mm heavy	1/2		75 violet				DS101PL		DS115PL	
30 mm	1/2		75 violet				DS18PL	DS19PL	DS20PL	DS72PL
30 mm heavy	1/2		75 violet				DS146PL		DS88PL	
35 mm	1/2		75 violet			DS27PL	DS28PL	DS29PL	DS30PL	DS31PL
			90 violet			DS22PL	DS23PL	DS24PL	DS25PL	DS26PL
35 mm heavy	1/2		75 violet					DS132PL	DS80PL	DS62PL
			90 violet				DS134PL		DS138PL	
37 mm heavy	1/2		90 violet					DS90PL	DS91PL	DS92PL

notes

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



**round needle with a cutting
edge (sharpened)**

needles		thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2	
length	curvature									
37 mm	1/2		75 violet					DS60PL	DS57PL	DS83PL
		75 undyed			DS119PL					
		90 violet			DS58PL	DS128PL	DS84PL	DS127PL	DS102PL	
		150 violet						DS162PL		
40 mm	1/2		75 violet				DS35PL	DS36PL	DS37PL	DS103PL
		90 violet				DS32PL	DS33PL	DS34PL	DS59PL	
40 mm heavy	1/2		75 violet					DS75PL	DS69PL	
		90 violet					DS133PL	DS86PL	DS79PL	
45 mm	1/2		75 violet				DS38PL	DS39PL	DS40PL	DS41PL
		90 violet						DS108PL		
45 mm heavy	1/2		75 violet				DS74PL	DS76PL	DS77PL	DS78PL
48 mm	1/2		75 violet					DS42PL	DS43PL	DS44PL
		90 violet					DS157PL		DS118PL	
48 mm heavy	1/2		75 violet							DS109PL
		90 violet						DS126PL		
48 mm	1/2		150 violet						DS124PL	DS142PL
80 mm	1/2		75 violet					DS53PL	DS107PL	
80 mm heavy	1/2		75 violet					DS130PL		

Catalogue contains the basic kinds of the connections of the needle with the thread. We can create non-standards connections acc. to individual orders.



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days, multifilament, coated



round needle with a cutting edge (sharpened)

needles		thread length (cm) color							
length	curvature		5/0	4/0	3/0	2/0	0	1	2
30 mm	hook	 75 violet						DS21PL	
31 mm	hook	 75 violet					DS148PL	DS98PL	
31 mm heavy	hook	 75 violet						DS120PL	

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round needle with a trocar edge

needles		thread length (cm) color								
length	curvature		5/0	4/0	3/0	2/0	0	1	2	
22 mm	1/2		75 violet			TS01PL	TS02PL			
			90 violet			TS03PL	TS04PL			
25 mm	1/2		75 violet			TS05PL	TS06PL	TS07PL	TS08PL	
			90 violet				TS09PL	TS10PL	TS11PL	
30 mm	1/2		75 violet					TS12PL	TS13PL	
			90 violet					TS14PL	TS15PL	
36 mm	1/2		75 violet					TS16PL	TS17PL	TS18PL
			90 violet					TS19PL	TS20PL	TS21PL
40 mm	1/2		75 violet					TS22PL	TS23PL	TS24PL
			90 violet					TS25PL	TS26PL	TS27PL

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round cutting needle (sharpened)

needles		thread	length (cm)	5/0	4/0	3/0	2/0	0	1	2
length	curvature									
22 mm	1/2		75 violet				RD02PL			
26 mm heavy	1/2		75 violet				RD01PL			
30 mm	1/2		75 violet						RD07PL	
35 mm	1/2		75 violet				RD06PL			
35 mm heavy	1/2		75 violet						RD04PL	
40 mm heavy	1/2		75 violet					RD03PL		
48 mm heavy	1/2		75 violet						RD05PL	

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



round separating needle

needles			thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
length	curvature									
17 mm	1/2		75 violet	RR09PL	RR10PL	RR08PL				
20 mm	1/2		75 violet		RR07PL	RR11PL				
22 mm	1/2		75 violet		RR06PL	RR03PL	RR04PL			
			90 violet				RR05PL			
26 mm	1/2		75 violet		RR12PL	RR02PL	RR01PL			
30 mm	1/2		75 violet				RR13PL			


notes



reverse cutting P needle blade

needles			thread length (cm) color	5/0	4/0	3/0	2/0	0	1	2
length	curvature									
19 mm	3/8		45 undyed		DM02PL	DM01PL				
22 mm	3/8		75 violet		DM04PL					
24 mm	3/8		75 violet				DM03PL			

notes



PGLA LACTIC (poly(glycolide-co-L-lactide))


surgical sutures, absorbable from 60-70 days,
multifilament, coated



reverse cutting needle and round needle

needles		thread							
length	curvature		5/0	4/0	3/0	2/0	0	1	2
30 mm	3/8	150 violet					DKORS01PL		
40 mm	1/2								
37 mm	1/2	150 violet					DKORS02PL		
37 mm	1/2								

notes



cutting needle and round needle

needles		thread							
length	curvature		5/0	4/0	3/0	2/0	0	1	2
30 mm	3/8	150 violet					DKRS01PL		
40 mm	3/8								

notes

PGLA LACTIC (poly(glycolide-co-L-lactide))

surgical sutures, absorbable from 60-70 days,
multifilament, coated

LIGATURES

threads length	color thread	7/0	6/0	5/0	4/0	3/0	2/0	0	1	2
2 x 70 cm	violet				P111PL	P114PL	P88PL	P66PL	P67PL	P89PL
2 x 75 cm	violet					P45PL	P46PL	P47PL	P48PL	P75PL
3 x 45 cm	violet			P90PL	P50PL	P01PL	P02PL	P03PL	P04PL	P53PL
3 x 75 cm	violet					P15PL	P16PL	P17PL	P18PL	P19PL
5 x 70 cm	violet					P68PL	P56PL	P57PL	P58PL	
6 x 45 cm	violet				P49PL	P05PL	P06PL	P07PL	P08PL	P09PL
6 x 75 cm	violet					P20PL	P21PL	P22PL	P23PL	P24PL
6 x 75 cm	undyed					P74PL				
7 x 75 cm	violet				P92PL	P51PL	P85PL	P91PL		
10 x 45 cm	violet				P65PL	P64PL	P63PL	P62PL	P61PL	
12 x 45 cm	violet				P44PL	P10PL	P11PL	P12PL	P13PL	P14PL
12 x 45 cm	undyed					P73PL				
75 cm	violet				P101PL	P60PL	P59PL	P82PL	P83PL	P84PL
120 cm	violet				P81PL	P76PL	P77PL	P78PL	P79PL	P80PL
140 cm	violet				P98PL	P99PL	P100PL	P93PL	P94PL	P95PL
150 cm	violet			P29PL	P30PL	P31PL	P32PL	P33PL	P34PL	P35PL
150 cm in reel	violet				Pr30PL	Pr31PL	Pr32PL	Pr33PL	Pr34PL	Pr35PL
250 cm	violet				P96PL	P36PL	P37PL	P38PL	P39PL	P40PL
250 cm in reel	violet				Pr96PL	Pr36PL	Pr37PL	Pr38PL	Pr39PL	Pr40PL
300 cm	violet						P113PL	P115PL	P116PL	P117PL
300 cm in reel	violet					Pr112PL		Pr115PL	Pr116PL	Pr117PL

notes